

Venus 3V & 4V

FEATURES

Disassemble in 1,2,3,4...5 parts



blue



6
km/h

Adjustable steering

Basket

optional
20A batteries

Swivelchair
adjustable in height

optional
20A batteries
available!

Full tires

VENUS 3V - 4V

Venus 4 sport



champagne



6
km/h

Frontlight

Basket

A compact 4-wheel scooter for indoor and outdoor usage. Comfortable and practical in day to day usage and can be disassembled in 5 parts without extra tools.



Full tires


Swivelchair
adjustable
in height





km/h



km



V3V	6	10,5	8°/14%	60	1100
V4V	6	10,5	8°/14%	60	1400
V4s	6	16	6°/10%	50	2470



V3V	510	900	1090	200 x 50 / 3	525-575	43	130	2 x 12V / 12 Ah	270 W
V4V	510	900	1070	200 x 50 / 4	525-575	46	130	2 x 12V / 12 Ah	270 W
V4S	550	900	1040	220 x 80 PU/4	500-550	55	115	2 x 12V / 20 Ah	270 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

