

VERMEIREN



N.V. Vermeiren N.V.

Vermeirenplein 1 / 15
B-2920 Kalmthout
Tel: +32(0)3 620 20 20
Fax: +32(0)3 666 48 94
website: www.vermeiren.be
e-mail: info@vermeiren.be

Vermeiren Austria GmbH

Schärdinger Straße 4
A-4061 Pasching
Tel: +43(0)7229 64900
Fax: +43(0)7229 64900-90
website: www.vermeiren.at
e-mail: info@vermeiren.at

Vermeiren Polska Sp. z o.o

ul. Łączna 1
PL-55-100 Trzebnica
Tel: +48(0)71 387 42 00
Fax: +48(0)71 387 05 74
website: www.vermeiren.pl
e-mail: info@vermeiren.pl

Vermeiren ČR S.R.O.

Nádražní 132
702 00 Ostrava 1
Tel: +420 596 133 923
Fax: +420 596 133 277
website: www.vermeiren.cz
e-mail: info@vermeiren.cz

Vermeiren Deutschland GmbH

Wahlerstraße 12 a
D-40472 Düsseldorf
Tel: +49(0)211 94 27 90
Fax: +49(0)211 65 36 00
website: www.vermeiren.de
e-mail: info@vermeiren.de

Vermeiren Italia S.R.L.

Viale delle Industrie 5
I-20020 Arese MI
Tel: +39 02 99 77 07
Fax: +39 02 93 58 56 17
website: www.vermeiren.it
e-mail: info@vermeiren.it

Vermeiren India Rehab Pvt. Ltd.

#400, Red Sanders Boulevard,
Sri City, Mallavaripalem (V), Sathyavedu
(M) Chittoor Dist., Pin - 517588
Tel: +91 80 41518899
Fax: +91 80 41528899
website: www.vermeiren.in
e-mail: info@vermeiren.in

Vermeiren Iberica, S.L.

Carretera de Cartellà, Km. 0,5
Sant Gregori Parc Industrial - Edifici A
17150 Sant Gregori (Girona)
Tel: +34 972 428 433
Fax: +34 972 405 054
website: www.vermeiren.es
e-mail: info@vermeiren.es

Vermeiren France S.A.

Rue de L'Empire
F-59133 Phalempin
Tel: +33(0)3 28 55 07 98
Fax: +33(0)3 20 90 28 89
website: www.vermeiren.fr
e-mail: info@vermeiren.fr

Vermeiren Suisse S.A.

Eisenbahnstrasse 62
3645 Gwatt (Thun)
Tel: +41 (0)33 335 14 75
Fax: +41 (0)33 335 14 67
website: www.vermeiren.ch
e-mail: info@vermeiren.ch

Vermeiren Medical Equipm. Co. Ltd.

18 Jiepu Road, SIP
215126 Suzhou
Jiangsu Province
China
website: www.vermeiren.com.cn
e-mail: info@vermeiren.cn

Vermeiren Nederland B.V.

p/a Vermeirenplein 1 / 15
B-2920 Kalmthout
Tel: +32(0)3 620 20 20
Fax: +32(0)3 666 48 94
website: www.vermeiren.co.nl
e-mail: info@vermeiren.co.nl



CONTENTS



05

MANUAL WHEELCHAIRS

43

ELECTRONIC WHEELCHAIRS

67

SCOOTERS - RAMPS

87

COMFORT CHAIRS

BEDS AND HOMECARE

99

BARIATRIC PRODUCTS

115

TRICYCLES

123

EASY LIFE RANGE

131



manual

Manual wheelchairs

Manual wheelchairs kids

708 kids.....	06
Eclips X4 kids.....	07
Sagitta kids	09
Eclips X4 90° kids	10
Gemini 2	11

Standard and light manual wheelchairs

708D - 708D T30	12
V100 - V200D - V300 - V300 30°	13
D100 - D200 - D200 30° - D400	17
Eclips+ - Eclips+ 30°	22
Eclips X4 - Eclips AD	24
Eclips X4 30° - 90°	26

Jazz S50	27
----------------	----

Manual active light wheelchairs

Sagitta	28
V300 XR.....	30
V500 - V500 30°	32

Stabilo armpads	34
V-drive	35
Options manual wheelchairs	36
Neoflex belts	38
Gel cushions.....	40
Bobby	42

708 kids

MANUAL KIDS WHEELCHAIR

Model 708 kids

Compact sturdy wheelchair in steel. Offset of the rear wheels for improved stability. Equipped with the Quickrelease system for fast and easy removal of the rear wheels.

FEATURES

Sturdiness on all levels



C16
dotted
blue



C29
grey

Supple nylon seat
and backrest

Soft padded
armrests

Quick-release
rear wheels

One-push button
system to detach
the armrests


Efficient
brake system

Tipping pieces
at the rear

Solid and durable
metal footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



500
550

310
310

320
370

210
210

445
445

420
420

210-310
210-310

350
350

985
985

18
18

60
60

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips X4 kids

MANUAL KIDS WHEELCHAIR

Ultralight kids wheelchair

Constructed in aluminum, the Eclips X4 kids offers you an ultra light wheelchair combined with sturdiness and ease of use. Equipped with clothing-protectors results in a very modern look. Without the wheels and the footrests, the Eclips X4 kids weighs no more than 9 kilograms making it very easy to transport.



model T30

Kit T30 push version with AR 12" wheels



Optional spoke protectors

Quick-release rear wheels

Aluminum front forks adjustable in height

Padded seat and back



Armrests to protect the clothes




Height adjustable handles (optional)

Internal and external collapsible footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



500
550

310
310

320
370

150-180
150-180

420-470
420-470

350-400
350-400

310/330/350/370
330/350/370/390

830-1000
830-1000

15
15

80
80

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Sagitta kids

MANUAL KIDS WHEELCHAIR

Active and ultralight

Active wheelchair for kids. The Sagitta kids is a mirror image of the Sagitta for adults. And even more important, this wheelchair grows with your child thanks to the variety of adjustments. The frame is constructed in such a way that you can change the seat width.




Adjustment of the push bar



Sport brakes (scissor brakes)

Wide range of accessories and options.

Refer to the accessories page (page 29 & 36 - 37) and check the compatibility table on page 62 - 65.



460-555	240	6°/11°	370-428 (with cushion)	3"	230-430	220-270	280/300/320	700-750	4°	9	60
480-575	260	6°/11°	410-468	4"	230-430	220-270	280/300/320	700-750	6°	9	60
500-595	280	6°/11°	(without cushion)	4"	230-430	220-270	280/300/320	700-750	8°	9	60
520-615	300	6°/11°			230-430	220-270	280/300/320	700-750	10°	9	60
540-635	320	6°/11°			230-430	220-270	280/300/320	700-750		9	60
560-655	340	6°/11°			230-430	220-270	280/300/320	700-750		9	60

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips X4 90° kids

LIGHT KIDS WHEELCHAIR

Light wheelchair for kids with a 90° backrest inclination. From the very first design, the Eclips X4 90° was made to offer a light, well designed wheelchair which is easy to handle and easy to transport.



C34
yellow

C35
pink

C36
lila

C37
blue



90° inclination



Optional B14 - Backrest extension



Padded seat and back


Quick-release rear wheels

Aluminum front forks adjustable in height

Internal and external collapsible footrests

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



500
550

310
310

320
370

150-180
150-180

420-470
420-470

350-400
350-400

310-330-350-370
330-350-370-390

830-1000
830-1000

16
16

60
60

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Gemini 2


KIDS BUGGY

The Gemini is an aluminum Reha Buggy and its design permits various possibilities of adaptation.


In this standard version we offer a:

- 5-points-security belt
- Spreader
- Head support

A buggy which parents adore not only for its excellent functionality, but also for the quality design and above all - the total weight is no higher than 13,5 kg.



Even for larger children we offer a model which is particularly designed towards stability. Manufactured from a steel frame construction to provide sturdiness and quality.



620
690

320
400

200
200

510
510

520-640
650-780

250-340
315-400

300-360
380-430

950
910

18
27

35
50

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

708D

STANDARD & LIGHT MANUAL WHEELCHAIRS

The armrests are foldable and removable and equipped with soft, extra long and comfortable padding. The fabric of the seat and back, together with aluminum legrests, will provide optimal comfort while using this chair.



C16
dotted
blue



C27
blue



C29
grey

FEATURES

Classic modular wheelchair



T30 version



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



560
590
610
630
650
670



280




390
420
440
460
480
500



225



490



420



383-458



430



1109



18
18
19
19
19
19



120

V100

FEATURES

Wide variety of possibilities

STANDARD & LIGHT MANUAL WHEELCHAIRS

The V100 is a light weight manually propelled wheelchair for temporary or permanent use with a range of settings.



C55
inox



C58
blue

- 3 seat heights
- 2 seat depths
- 3 armrest heights
- 3 depths for the arm pads
- wide variety of options

- Ergonomical push handles
- Height adjustable


Foldable and detachable armrests

Foldable and detachable footrests

Rear wheels with a Quickrelease system

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



570		390								
600		420								
620	300	440	220-240	470/500/530	420	345-515	430/450	1010	18	130
640		460								
660		480								
680		500								

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V200 D

STANDARD & LIGHT MANUAL WHEELCHAIRS

The V200 D is a light weight manually propelled wheelchair for temporary or permanent use with optimised settings.



C59
champagne



C60
brown

- 4 seat heights
- 2 seat depths
- 5 armrest heights
- 5 positions for the armrests
- Mobile phone pouch under the seat
- A range of available accessories

FEATURES

A multitude of adjustments and a wide variety of options



Angelwings backrest, providing comfort for specific needs for your back

- Pushhandles adjustable in height
- String bag

In 3 heights and depths adjustable armrests

Foldable and detachable footrests

Rear wheels with Quickrelease system

Angle adjustable footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



590		390								
620		420								
640	320	440								
660		460								
680		480								
700		500								
			220-240	440/470/500/530	400	345-510	440/460	1010	16,5	130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V300

FEATURES

A multitude of adjustments
with a seatwidth up to 600 mm

STANDARD & LIGHT MANUAL WHEELCHAIRS

The V300 is a light weight manually propelled wheelchair for temporary or permanent use with maximized settings.



C59
champagne



C80
blue

Adjustable seat depth
by shifting the back



Tension adjustable backrest



Adjustable seat depth

In height and depth
adjustable armrests

Tension adjustable
back and seat

Foldable and
detachable armrests

Foldable and
detachable
footrests

Rear wheels with
Quickrelease
system

Angle
adjustable footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



590	390									
620	420									
640	440	220-240	440/470/500/530	420	850/880/910/940	440/460/480/500/520	1040	17,7	130	
660	460									
680	480									
700	500									
740	530				840/930	490/510/530	1085	21	170	
770	560									
810	600									

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V300 30°

STANDARD & LIGHT MANUAL WHEELCHAIRS

- 4 seat heights
- 5 seat depths
- 3 armrest heights
- 6 seat widths
- Armpads adjustable in depth
- Mobile phone pouch under the seat
- A range of available accessories



C59

champagne

Removable pushbar

FEATURES

A multitude of adjustments
with a 30° reclinable backrest



Optional preshaped back and seat L34 - L35

30° backrest
inclination

Special designed armrests which are
adjustable, foldable and removable

Angle
adjustable footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



590	390									
620	420									
640	440	220-240	440/470/500/530	420	850/880/910/940	440/460/480/500/520	1040	17,7	130	
660	460									
680	480									
700	500									

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

D100

FEATURES

Basic comfort and light to transport

STANDARD & LIGHT MANUAL WHEELCHAIRS

The D100 is your ideal compact lightweight wheelchair. Completely foldable and plenty of adjustment possibilities towards your own desires. Optional **BZ8** available.



C70
grey

Pushhandles
adjustable in height

Seat height
adjustable in
3 positions

- Foldable and
detachable armrests

- In height and depth
adjustable armpads

Foldable and
detachable
footrests

Integrated lowering
of the rear axle

Seat height adjustable
in 3 positions

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



580	380									
600	400									
620	420									
640	440	220-240	470/500/530	420	880/910/940	430	1010	18	130	
660	460									
680	480									
700	500									
720	520									

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

D200

STANDARD & LIGHT MANUAL WHEELCHAIRS

The D200 is your ideal hybrid lightweight aluminum wheelchair. Completely foldable and lots of adjustment possibilities, including the footplates which can be adjusted in angle.



C70
grey

Pushhandles
adjustable in
height

Seat height
adjustable in
3 positions

- Foldable and
detachable armrests

- In height and depth
adjustable armpads

Foldable and
detachable
footrests

Integrated lowering
of the rear axle

Angle
adjustable footplates

FEATURES


Easy, foldable and light



*Model D200 and D100 are also available
with an optional B69 (foldable) backrest*

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



580
600
620
640
660
680
700
720

300

380
400
420
440
460
480
500
520

220-240

470/500/530

420

385-455

430

1010

16,8

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

D200 30°

FEATURES

Light and functional
with a 30° reclinable backrest

D200 30°

STANDARD & LIGHT MANUAL WHEELCHAIRS

The D200 30° is a lightweight aluminum wheelchair with 30° reclining backrest and foldable to the front. It also offers many possibilities and adjustments.



C70

grey

- 30° reclinable backrest
- Metal push bar




Optional BZ7 - BZ8 legrests



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



580		380								
600		400								
620		420								
640	300	440	220-240	470/500/530	420	385-455	430	1010	16,8	130
660		460								
680		480								
700		500								
720		520								

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

D400

STANDARD & LIGHT MANUAL WHEELCHAIRS

The D400 is your ideal compact lightweight wheelchair which is foldable with numerous adjustment possibilities. Tension adjustable back, different seat heights and seat depths... the D400 allows itself to be adjusted to your needs.



C70
grey

Pushhandles
adjustable in height

Tension adjustable
back



Foldable and detachable armrests
In height and depth adjustable armrest pads

Individual seat position:

- Seat heights
- Seat depths

Foldable and
detachable
footrests




Fork adjustable
in angle

Seat height adjustable
in 3 positions

Angle
adjustable footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



580
600
620
640
660
680
700
720

310

380
400
420
440
460
480
500
520

205-280

440/470/500/530

400-460

330-500

440/460/480/500/520

1040

17,7

130



Eclips+

STANDARD & LIGHT MANUAL WHEELCHAIRS

Lightweight aluminum wheelchair with a fixed back and reinforced cross. Equipped with Quickrelease rear wheels, foldable and removable armrests.



C25
diamond
black



C29
grey



C80
blue

FEATURES

A solid classic throughout the years



Optional (HEM2) hemiplegic kit



Optional (T30) push version



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



590
620
640
660
680
700
750



300




390
420
440
460
480
500
550



220



430-530



400-460



390-460



420-440-460



1040-1080



kg

14,5
14,5
15
15,8
15,8
15,9
16



kg

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips+ 30°

STANDARD & LIGHT MANUAL WHEELCHAIRS

Lightweight aluminum wheelchair with a reclining backrest of 30°, reinforced cross and Quickrelease wheels. Standard equipped with short armrests for easy table access.



C30
grey

FEATURES

Swift 30° adjustment of the backrest



Height adjustable handles



Inclination of the back by hydraulic springs for more comfort during the tilting



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Eclips X4

STANDARD & LIGHT MANUAL WHEELCHAIRS

Light aluminum wheelchair with adjustable frame. Adjustable in depth, reinforced cross and Quickrelease rear wheels. The footrests can be folded both inward and outward. For optimum comfort, the footrests are adjustable in depth.



FEATURES

Longlasting classic with a multitude of adjustments



Also available as 30° and 90° model (see page 26)



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



590 620 640 660 680 700 740	300	390 420 440 460 480 500 550	220	390-520	420-460	345-460	420-500	1000-1210	17 17,1 17,1 17,1 17,2 17,4	150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips AD

STANDARD & LIGHT MANUAL WHEELCHAIRS

The Eclips AD is equipped with a telescopic cross that allows adjusting the width of the seat so that one wheelchair can cover seat widths from 390 to 500 mm.



C29
grey

This light wheelchair comes with a fixed comfortable nylon backrest and aluminium, height adjustable, 24" rear wheels with a Quickrelease fixation. The front wheels, 200 x 35 mm, are made out of Krypton. Additionally the legrests can be removed or folded in- or outwards and can be adjusted in depth. As for the armrests, they can fold backwards or be totally removed for easy transfers.

FEATURES

6 different widths
combined in one wheelchair



Regulate the width of the base



590		390								
620		420								
640	300	440	220	430-530	420-460	390-460	420	1040	17	130
660		460								
680		480								
700		500								

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips X4 30° / 90°

STANDARD & LIGHT MANUAL WHEELCHAIRS

Eclips X4 30° or 90° - lightweight aluminum wheelchair with a reclinable backrest of 30° or 90°.
Ideal for longer term use of a wheelchair with a wide variety of options to fit the needs of the user.



C30
grey



C80
blue

Removable metal
push bar

Hydraulic springs

Optional B05 armrests - Foldable,
removable and height adjustable

Optional foldable and removable
footrests with calf support,
adjustable in angle and length

Rear wheels with
Quickrelease system

Foldable and
length adjustable
footplates



Eclips X4 30° with nylon seat and backrest

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



590
620
640
660
680
700
740



300




390
420
440
460
480
500
550



220



390/445/471/520 420-460




900-1015 420/440/460/480/500



1000-1210



19,5



130

Jazz S50

FEATURES

Your basic foldable wheelchair

JAZZ S50

STANDARD & LIGHT MANUAL WHEELCHAIRS

The wheelchair Jazz S50 is an excellent combination between technicality, comfort and Vermeiren quality.



C27
blue



--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Sagitta

MANUAL ACTIVE LIGHT WHEELCHAIR



C20
blue



C41
green



C52
gold



C53
orange



C55
inox



C56
blue



C76
white



Limited size
for transport

Option
inclinable backrest 6° - 14°

Seat and back made
from strong durable nylon

Tension-adjustable back and seat

Adjustable wheel angles:
0°, 2°, 4°, 6°, 8°
(order only)



Push handles
adjustable in height



Option B23 - Extra cushion



Option B78
Pair of anti-tipping wheels



Option B82 - Clothing protectors

Wide range of accessories and options.

Refer to the accessories page (page 29, 36 - 37) and check the compatibility table on page 62 - 65.



550-680
570-700
590-720
610-740
630-760
650-780



360
380
400
420
440
460




7°
7°
7°
7°
7°
7°



480
(with cushion)
510
(without cushion)



3"
4"
5"



230-430
230-430
230-430
230-430
230-430
230-430



100
adjustable



330-380-420
(standard)
440-460
(long version)



755-855
755-855
755-855
755-855
755-855



0°
2°
4°
6°
8°



>9,5




100

Options



Clothing protectors B82S



Clothing protectors B82



Clothing protectors
B82S + tube



Sport brakes
(scissor brakes)



Standard brakes



Footrest with
aluminium plate (kids)



Footrest with
separated footplates



Monoblock footrest
standard tube




Monoblock footrest
with aluminum plate



Push handles
adjustable in height



Standard
push handles



Wheels 3" translucent



Wheels 4" translucent



Wheels 5" translucent



Wheels 3" black



Wheels 4" black



Wheels 5" black



Spoke protector - B85
M3 "Extreme"



Spoke protector - B85
M1 "Flowers"



Spoke protector - B85
M2 "Flowers"



Spoke protector - B85
M3 "Pirate"



Spoke protectors - B85
logo Sagitta



Spoke protectors - B85
(transparent)



Delphi Alu Schwalbe Right Run

About 1,38 kg with tire
Compatible with tires:
1. Marathon Plus Evolution (1901721)
3. Schwalbe Right Run (2990948)



24" Wheels Spinergy

About 1,79 kg with tire
Compatible with tires:
1. Marathon Plus Evolution (1901721)



24" Wheels Off Road

About 2,24 kg with tire
Compatible with tires:
2. Landcruiser (2991487)



24" Wheels Viking

About 1,85 kg with tire
Compatible with tires:
1. Marathon Plus Evolution (1901721)

V300 XR

MANUAL ACTIVE LIGHT WHEELCHAIRS

FEATURES
Striking colours and design



C71
blue

C72
red

C73
green

C80
blue



Footplates adjustable
in angle



wheelchairs displayed with optional Viking rear wheels.

Options



Wheels 22"x 1 3/8" alu + wheel
rims in anodised aluminium
Schwalbe Right Run



Wheels 24"x 1 1/8" alu
+ wheel rims in anodised
aluminium
Schwalbe Right Run/DELFI



Wheels 24" Viking



Wheels 24" Spinergy Schwalbe
Marathon Plus Evolution

- 4 seat heights
- Easy and fast removal
of the rear wheels



Wheels 4" black



Wheels 5" black



Push handles
(adjustable in height)



Anti-tipping wheels



Travel wheels



Monoblock footrest
with aluminum plate



Manual articulating
legrest

V300 XR

- 4 seat heights
- 5 seat depths
- 9 seat widths (from 390 mm to 600 mm)
- 5 seat inclinations from 0° to 10°
- Mobile phone pouch under the seat
- Wide variety of options and features

Adjustable seat depth
by moving the backrest
and the cross



Seat and backrest tension adjustable



Internal and external
collapsible footplates

Wide range of accessories and options.

Refer to the accessories page (page 30, 36 - 37)
and check the compatibility table on page 62 - 65.



390
420
440
460
480
500
530
560
600

440-470-500-530

300
350
400

4"

5"

830-930

440-460-480-500-520
440-460-480-500-520
440-460-480-500-520
440-460-480-500-520
440-460-480-500-520
440-460-480-500-520
490-510-530
490-510-530
490-510-530

1040
1040
1040
1040
1040
1040
1085
1085
1085

> 14*

130
130
130
130
130
130
170
170
170

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V500

FEATURES

Light and trendy

STANDARD & LIGHT MANUAL WHEELCHAIRS

The V500 is a manual lightweight wheelchair for temporary or permanent use with maximum adjustments.



C80
blue



C81
red

Tension adjustable
backrest and seat

Design and lightweight
removable armrests.
Adjustable in height
and in depth

Folding brake
handles

Adjustable
seat depth

Removable and
foldable footrest
(in and outwards)

Anti-tipping wheels

Rear wheels with
Quickrelease system

Adjustable
seat depth

Angle
adjustable footplates

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



590
620
640
660
680
700

390
420
440
460
480
500

440/470/500/530

400

150x30

850/880/910/940

440/460/480/500

1010

15,6

130

V500 30°


FEATURES

V500 with easy
reclinable backrest

V500 30°

STANDARD & LIGHT MANUAL WHEELCHAIRS

The V500 30° wheelchair is a lightweight manual wheelchair with all the highlights of the standard V500 (modern design, seat and backrest tension adjustable, brakes with foldable and long handles, anti-tipping wheels, wide variety of adjustments...) with the added benefit of a 30° reclining backrest with hydraulic springs for optimum precision and comfort for the user as well as for the accompanying person.



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



590
620
640
660
680
700

390
420
440
460
480
500

440/470/500/530

400

150x30

850/880/910/940

440/460/480/500

1010

15,6

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Stabilo comfort armpads

Hemiplegic cuff allowing to form and reform the cuff as many times as you want and thus to give it the desired shape.

Ideal for patients having deformations of the lower arm.
This shapable armpad adapts itself perfectly to the desired shape, releasing tension and pressure from the hand or the wrist.

cm	width	width	T*
M	25	30	4
L	25	45	7

*T = thickness

Adjust according to your needs ...



Available
in two sizes



V-drive

Comfort

The V-Drive system simplifies your movements with a manual wheelchair. The V-Drive takes over from your pushing assistant. Taking a steep hill or inclinations, road bumps or just a small journey outside... all this is ideal for the V-Drive.

Standard version
for people up to 135 kg
(weight chair included)

H.D. version
for people up to 200 kg
(weight chair included)



Assembly

A very simple installation for the left or right side mounting of the V-Drive. Disassembly of your V-Drive for transport inside your car can be done in a minute.

Disconnect the battery pack and the driving system and you are ready to go.



New height adjustment to keep the motor horizontal

Delivered in a transport bag



standard
H.D.

390-490
500-600

6 km/h
6 km/h

8°
8°

15 km
15 km

2 x 12 V 12 Ah
2 x 12 V 12 Ah

150 W
200 W

21
17,4

135
200

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Options

MANUAL WHEELCHAIRS

HEADRESTS



B14 - backrest extension



B47 - adjustable headrest



L55 - shapeable headrest

SEATS & BACKRESTS



B17 - toilet seat



B22 - spreader



B23 - cushion 40 mm



B26 - cushion 20 to 60 mm



B37 - 60 mm cushion with oval opening



B42 - shopping net



B73 - arthrodesis seat



B79 - document holder




B90 - push bar



L04 - side support



L14 - anatomic back (soft)




L15 - anatomic seat (soft)



L14 - L15 - anatomic back (skai)



L15 - anatomic seat (skai)



L24 - preshaped back (skai)



L25 - preshaped seat (skai)



L34 - preshaped back (soft)



L34 - preshaped seat (soft)



L44 - tension adjustable backrest

ARMRESTS



B02 - long armrest



B03 - short armrest



B05 - armrest adjustable in height



B112 - hemiplegic armrest



B66/o - gutted rotating armrest



Bi5 - adjustable in height

KITS



AMP - kit: offset of the rear wheels



HEM2 hemiplegic kit: double rims



S1 - lowered seating kit



T30 - pushing wheelchair kit

Options

MANUAL WHEELCHAIRS



B06 - standard footrest



B18 - heel support



B19 - foot straps



B21 - wide calf support XL



BA6 - shortened footrest



BB6 - adjustable short footrest



BK6 - one-piece footrest



BL6 - one-piece footrest shortened



BN6 - shortened footrest



B07 adjustable footrest - long



B08 adjustable footrest - short



BK7 / BL7 AMP legrest



BP7 - plaster leg legrest



BZ7 / BZ8 - articulated legrest self-correcting in length
07 long - 08 short



BZ9 - one-piece adjustable legrest



BZ10 - BZ9 legrest with side corners of 70 mm



B12 - wooden table



B20 - safety belt



B31 - crutch holder



B33 - rim cover



B40 - brake extenders



B52 - serum holder



B58 - safety belt auto lock



B74 - brakes



B78 - anti-tipping



B83 - cannister holder



B85 - spoke protectors




B86 - travel wheels



B91 - anti-theft slot



B101 - anti-theft system with money-coin



B115 - plexi table

FOOTRESTS

LEGRESTS

ACCESSORIES

Neo-Flex

Range of Neo-Flex safety harnesses. 11 solutions completely compatible and adjustable to our wheelchairs. These harnesses combine comfort, solidity and reliability. Each model is available in several widths for the comfort of each user.



U73 - Harness with zip
Available in narrow or wide version
Available in 5 sizes
(XS, S, M, L and XL)



U74 - Spider harness without zip
Available in narrow or wide version
Available in 5 sizes
(XS, S, M, L and XL)




U76 - Abdominal belt - 2 attachment points. Available in 4 sizes (XS, S, M and L) and exists in 4 points of attachment version (**U75**)



U77 - Ankle support (pair)
Available in 5 sizes
(XS, S, M, L and XL)



U78 - Basic shoulder harness -
2 ways to attach it to the front (crossed or uncrossed)
Available in 3 sizes (S, M and L)



U79 - H-type shoulder harness
Available in 4 sizes (S, M, L and XL)



U80 - Harness waist retention
Available in 5 sizes (XS, S, M, L and XL)



U01 - Pelvic belt
Available in 4 sizes (XS, S, M and L)

Neo-Flex

For every problem Neoflex brings you the solution.
Other examples of use are available on our website www.vermeiren.com

Some
examples
of use

KYPHOSIS



LORDOSIS



SCOLIOSIS



ROTATION OF THE TRUNK



Gel cushions



V10 - protective cover
60 x 240 mm (pair)



V12 - V14 - Protective gel covers
fortubing 6" - 8" (pair)
120 x 160 mm - 120 x 200 mm



V16 - Gel heel protectors (pair)
120 x 140 mm



V18 - Velcro gel pad for walker
or crutch (pair) - 120 x 130 mm



V20 - Gel protection for a
safety belt (2 parts)
325 x 80 mm and 230 x 70 mm



V22 - Gel protection for gutted
armrest (pair) - 110 x 340 mm



V24 - Gel leg side protection
(pair) - 310 x 90 mm



V26 - Gel strap support for legs
2 pieces - 400 x 190 mm



V28 - Gel plate footrest
(pair) - 155 x 210 mm



V30 - Adhesive oval pad
(pair) - 80 x 100 mm



V32 - Adhesive oval pad
(pair) - 100 x 200 mm



V34 - Adhesive oval pad
(pair) - 120 x 300 mm



V36 - Adhesive ovoid pad
120 x 70 mm



V38 - Gel armrest pad
(pair) - 90 x 300 mm



V40 - Gel protective knee pad
(pair) - 100 x 100 mm



V10



V14




V12




V18



V22



V32



V28



V10

V12

Gel cushions

In many cases, there is a need for additional protection of certain body parts during the use of a wheelchair. Vermeiren is proud to offer a wide range of 16 types of comfort gel cushions that are very easy to install and can be used in almost any situation.



V24



V26



V38



V40



V30



V32



V34



V36

Bobby

TRANSPORT WHEELCHAIR

The Bobby is your ideal compact, foldable transferchair. Light and easy to fold in seconds. Either on holiday, hospital or airport for a transfer from point A to point B, the Bobby chair is easy to fold and to store and requires little maintenance.



C80

blue

- Foldable
- Two manual brakes
- Padded armrests
- Angle adjustable footplates
- Adjustable heelstraps




FEATURES

Compact - Sturdy - Functional



Option: BZ7- BZ8 articulated legrest self-correcting in length (07 long, 08 short)



590
640

270
270

420
480

220
220

500
500

380
380

420
420

950
950

9,40*
9,95*

115
115

* In optimal configuration without legrests. Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

electro

Wheelchairs

Electronic wheelchairs kids

Springer	44
Forest kids	45

Electronic wheelchairs adults

Timix	46
Timix S.U.	47
Navix (RWD)	48
Navix S.U.	49
Tracer	50

Express	51
Squod	53
Forest 3	54
Forest 3 S.U.	55
Forest 3 Initial	56
Forest 3 Plus	57
LiNX controller	58
Options	60

Springer

ELECTRONIC KIDS WHEELCHAIR

The Springer is a perfectly designed electronic wheelchair for kids allowing maximum adjustments for a complete personal touch. The Springer can be used indoor and outdoor on most terrains. Complete with LED lighting at the front and the back, adjustable armrests and backrest and adjustable leg- and footrests. Even further, you can adjust the seat width and seat depth for optimal personal comfort.



C67
orange

Anti-tipping
wheels
at the back



Profiled wheels for a
better grip on uneven terrains

Ground clearance 80 mm

LED lights in front and back

4 shock absorbers

Seat LIFT up to 300 mm for

- Better social interaction at similar level with friends & family at home, at school, ...

Better access to places difficult to reach


- Shelves at school, in library, at home, in telecom shop, ...
- Cupboard, drawers,... in kitchen, living room, bath & bed room, ...
- Easy outdoor & indoor manoeuvring
- Seat lift for better access & social interaction
- Great looking & highly configurable



Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



2x12V-50Ah	2x220W 2x350W	6 km/h 10 km/h	23 km	6"/10%	70	1400
------------	------------------	-------------------	-------	--------	----	------



600	260-380	160-270	530-810	310-370	370-490	15°	-5/+55°	1010	90	50
-----	---------	---------	---------	---------	---------	-----	---------	------	----	----

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Forest kids

FEATURES

Kids wheelchair
for outdoor terrain

ELECTRONIC KIDS WHEELCHAIR

The Forest Kids is an electronic wheelchair for children, used for indoor and outdoor and ideal in overcoming obstacles. The 4 wheels are always in contact with the ground even on rough surfaces up to 120 mm.



C53
orange



C55
inox



C57
red



C71
blue



Electronic tilting of seat and back



Height adjustable
pushing handles

Armrests adjustable
in height

Easy to operate
controller

Legrests adjustable
in depth

LED lights in front
and back



Legrests adjustable
in width




Stability axle



2x12V-75Ah 2x220W 2x350W 6 km/h 10km/h 30 km 10°/17% 60 1400

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



600 320 350 370 200-250 530 330-390 350 15° -5°/55° 1050 95 50

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

ELECTRONIC ADULT WHEELCHAIR

The TIMIX allows for sharper turning indoor and for better overcoming obstacles like thresholds when entering the local grocery store. The user friendly TIMIX is adjustable to individual needs and allows for easy transfers to bed, bath, toilet, car etc...



C71
blue



C72
red



C73
green

- Full suspension on all wheels
- Stability axle ensuring contact of all four wheels with the ground on uneven terrains up to 90 mm
- Front and rear lights and indicators



Optional 300 mm lift

- Kitchen cupboards
- Cash dispensers
- Bar
- etc...

- Adjustable seat depth
- Optional luggage rack

- Multi-functionally adjustable armrests, swing-away and removable
- Electrically adjustable back and seat
- In- and outwards swing-away controller
- Controller adjustable in height
- No pinch points: safety for fingers & clothes

Fully integrated motor release levers

Anti-tipping wheels



2x12V-80Ah	2x220W or 2x350W	6 km/h 10 km/h 14 km/h	63 km 40 km 39 km	9°	100	1320
------------	------------------------	------------------------------	-------------------------	----	-----	------

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



640	400	242-350	490/540	420-520	630	5°-15°	1270	125	130
640	450	242-350	490/540	420-520	630	5°-15°	1270	125	130
660	500	242-350	490/540	420-520	630	5°-15°	1270	125	130

Timix s.u.

ELECTRONIC STANDUP WHEELCHAIR

The Timix S.U., also known as the Timix Stand Up, is an electronic stand-up wheelchair completely controllable from its controller. From a vertical stand-up position to a complete horizontal position in just one movement.



C71
blue



C72
red



C73
green



Fully collapsible knee-lock



Chest safety belt
with velcro

Back and seat
in soft version

Automatic length
compensating footplate

Anti-tipping wheels



2x12V-80Ah	2x220W 2x350W	6 km/h 10 km/h 14 km/h	63 km 40 km 39 km	9°	100	1320
------------	------------------	------------------------------	-------------------------	----	-----	------

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



670

450

min. 242
&
max. 350

540
(without cushion)

450
500

610

1260

160

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Navix (RWD)

ELECTRONIC ADULT WHEELCHAIR

The Navix allows for sharper turning indoor and for better overcoming obstacles like thresholds when entering the local grocery store. The user friendly Navix is adjustable to individual needs and allows for easy transfers to bed, bath, toilet, car etc..



C30
metallic
grey



C71
blue



C72
red



C73
green

**REARWHEEL
DRIVE
available**

FEATURES
Better manoeuvring with
a compact front-wheel drive



- Better overcoming obstacles
- Better taking sharp turns
- Easy to reassemble
- Total width for Navix limited to 570 mm
- Ground clearance of 60 mm
- Adjustable to your needs

Better grips on outdoor surfaces
with profiled large front wheels



Optional electrical inclination
of the seat and backrest



Optional 300 mm lift for
a better accessibility



Optional: BZ7-E
electrical adjustable legrests



2x12V 38Ah 2x220W 6 km/h 20 km 9° 60 1250

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



580	400	260-350	400/450/500	380-500	580	330-400	1010 (945 RWD)	79 kg	130
580	450	260-350	400/450/500	380-500	580	330-400	1010 (945 RWD)	79 kg	130
630	500	260-350	400/450/500	380-500	580	330-400	1010 (945 RWD)	79 kg	130

RWD: Rear wheel drive - Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Navix s.U.

ELECTRONIC STANDUP WHEELCHAIR

The Navix S.U., also known as the Navix Stand Up, is an electronic stand-up wheelchair similar to the standard Navix with added functionality.



C30
metallic
grey



C71
blue



C72
red



C73
green



Fully collapsible knee-lock



2x12V-50Ah 2x220W 6 km/h 32km 9° 50 1250

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



670

450

min. 242
&
max. 350

470

450
500

610

1090

122

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Tracer

ELECTRONIC ADULT WHEELCHAIR

The Tracer is equipped with a control unit with microprocessor that gives an extra dimension in using an electronic wheelchair: it is easy to control, with a minimum of effort and maximum safety. The motors in the Tracer are the latest generation of motors. 2 x 200 Watts enables you to go up slight inclines without problems.



C20
blue



C21
red



C30
grey



C44
orange

One piece metal pushbar



Anti-tipping wheels

Anatomical seat and back with a safety belt

Front and rear light with indicators

FEATURES

Available with 550 mm seat width

model TRACER GT

- Compact outdoor wheelchair
- Electrical seat and back
- Nylon seat and back



Colour LCD display



Optional: seat and back in Soft version



2x12V-75Ah 2x200W 10 km/h 35 km 10°/17% 70 1400

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



660-720	390	200-240	550	450	600	360-510	1190	119	125
680-780	450		550	450	600		1190	119	125
770-870	550		580	500	600		1190	128	150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Express

ELECTRONIC ADULT WHEELCHAIR

Electronic wheelchair with a double cross, equipped with footrests and armrests adjustable in height and width. Complete with front wheels 200x50 and rear wheels 12¹/₂ x 2¹/₄. Equipped with a powerful engine of 2 x 200W and 38 A/h batteries.



Easy and user-friendly controller

FEATURES

Disassemble fast and easy

LINX[™] ACCESS

Wheelchair Control System

LE Controller System :

- Ideal for light weight, indoor wheelchairs
- Easy to use
- Simple
- Reliable



Very compact and foldable chair which facilitates transport



- Adjustable and removable footrests to the in- and outside.
- Footplates adjustable in angle



2 separated batteries



2x12V-38Ah	2x150W	6 km/h	18 km	7°	15	1400
2x12V-50Ah			25 km			

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



580		390								
610		420								
630	400	440	min. 225	540	430	430	400-520	950	67	90
650		460	max. 315							
690		500								

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

LiNXTM

The next generation wheelchair control system

- A seamless and expandable system
- Modular components – mix and match modules
- State-of-the-art technology
- Revolutionary approach to programming



Wireless Programming & Diagnostic Tools for PC & iOS Devices

- Dedicated access levels for manufacturers and their distributors.
- Bluetooth wireless communications to the PC or iOS devices.



Drive statistics Live diagnostics

All the information you need at your fingertips on your mobile devices*.



* LiNX Access iOS - Compatible with:

- iPhone 5 and newer
- iPod Touch Gen 5 and 6
- iPad 2 and newer including iPad Pro and all iPad mini's

Squod

ELECTRONIC ADULT WHEELCHAIR

The Squod is a perfect electric wheelchair for both inside and outside usage. The seat and backrest are made in Soft, which gives you the opportunity to feel more comfortable and at ease in any situation. The joystick can fold to the side making it more convenient if you want to drive closer to a table.

C30
diamond
grey

C41
green

The seat and backrest are
in Soft for a more
comfortable experience



Reclinable backrest
up to 30°



Anti-tipping wheels



Armrests can be adjusted
in height and width



Removable armrests
Adjustable in height
and width

- Adjustable and
removable footrests to
the in- and outside.
- Footplates adjustable
in length

SQUOD



2x12V-70Ah 2x200W 6 km/h 30 km 6°/10,5% 70 1400
2x12V-38Ah width 390

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



590 390-430 200-250 530 430 510 1070 111 130
630 440-490
700 500-550

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Forest 3

ELECTRONIC ADULT WHEELCHAIR

The FOREST 3 is ideal for crossing obstacles in front of you. The 4 wheels of the Forest will always be in contact with the ground in spite of uneven terrain up to 100 mm in height differences.



New LiNX controller system



Optional LIFT function available



C71

blue

C72

red

C73

green



C76

white

C77

black

C78

blue

Folding armrests, adjustable in height, in depth, in width and inclination with position memory

Ventilated coating for the seat and the back avoiding any problem of perspiration

Electrical adjustment of the seat and backrest

Seat height starting from 415 mm (with electrical seat adjustment)



Suspension on the rear wheels

Ground clearance of 75 mm

Stability axle in front



2x12V-80Ah

2x220W
2x350W

6 km/h
10 km/h

50km
42km

9°

100

1320

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



595
610
650

400
450
500

242-350
242-350
242-350

430-480
(sur commande)

420-520
275-395
420-520

630
630
630

275-395
1185
275-395

1185
1185
1185

126
126
126

130
130
130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Forest 3 s.u.

FEATURES

Multiple positions by the touch of a button

ELECTRONIC STANDUP WHEELCHAIR

The Forest 3 S.U., also known as the Forest 3 Stand Up, is an electronic stand-up wheelchair similar to the standard Forest 3 with added functionality.



C71
blue

C72
red

C73
green



C76
white



C77
black



C78
blue



Fully collapsible knee-lock

Swing away controller to left or right and adjustable in height

Upholstered seat and back

Automatic length compensating footplate

2x2 wheels for optimal support



2x12V-80Ah 2x220W 6 km/h 50 km 9° 100 1320

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



670

450

min. 242
&
max. 350

530

min 450
&
max 500

610

1170

165

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

FOREST 3 S.U.

Forest 3 initial

FEATURES
Forest 3 technology
in its purest form

ELECTRONIC ADULT WHEELCHAIR

Ideal for indoor and outdoor use, the Forest 3 Initial is a refined version of his big brother, the Forest 3 but keeps the main strengths facilitating in particular the crossing of obstacles.

The front stability axle allows the Forest 3 Initial to keep his 4 wheels on the ground despite surface differences up to 10 cm. Presented in a 6 km/h version and equipped with the new **LiNX LE** control system.

LiNX™ ACCESS
Wheelchair Control System


LE Controller System :

- Ideal for light weight, indoor wheelchairs
- Easy to use
- Simple
- Reliable



2x12V-80Ah 2x220W 6 km/h 50km 9° 100 1320

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



595	400	242-350	430-480	420-520	560	275-395	1185	133,5	130
610	450	242-350	(on order)	420-520	560	1185	1185	133,5	130
650	500	242-350		420-520	560	275-395	1185	133,5	130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Forest 3 Plus

ELECTRONIC ADULT WHEELCHAIR

The Forest 3 Plus is a wheelchair that covers most of the features of a normal Forest 3 but specially designed for larger sizes.



C71
blue



C72
red



C73
green



C76
white



C77
black



C78
blue

250KG

FEATURES

A special Forest 3 for wider sizes



New LINX controller system

LINX™ ACCESS
Wheelchair Control System

LE Controller System :

- Ideal for light weight, indoor wheelchairs
- Easy to use
 - Simple
 - Reliable



2x12V-80Ah 2x350W 6 km/h 49 km 9° 52 1840

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



705
760

500-550
600-650

335-380
335-380

450
450

450/500/550
450/500/550

590
590

250-399
250-399

1145-1163
1145-1163

134,35
134,35

250
250

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



THE CONTROLLERS FOR A NEW GENERATION OF WHEELCHAIRS

LiNX™ ACCESS
Wheelchair Control System



LiNX LE controller with basic functions

LiNX™



Standard LiNX controller with many possible functions such as inclination of the seat, backrest, lights etc ...


Standard controller on the wheelchair

Compatible with LiNX LCD

FOREST 3	LiNX	yes
FOREST 3 +	LiNX	yes
FOREST 3 Initial	LiNX LE	no
FOREST 3 S.U.	LiNX	yes
NAVIX FWD	LiNX	yes
NAVIX RWD	LiNX	yes
NAVIX S.U.	LiNX	yes
TIMIX	LiNX	yes
EXPRESS	LiNX LE	no
SPRINGER	LiNX	yes
FOREST kids	LiNX	yes

LiNX™

with **LCD touch screen**



NEW

LiNX control with LCD touch screen.
Allows access to a large number of functions (inclination of seat and backrest, lighting ...) with direct visualization on the LCD screen for more comfort and practicality.

Simple and robust

- Intuitive control
- Ergonomic joystick with minimal force required (190g)
- Simultaneous driving and seat adjustment
- Wireless connection with smartphone or tablet

Adapts to different types of terrain

LiNX steers very responsive on any terrain and even at low speed. The controller will continuously follow up on and adapt to the underlay. This enhances the control you have over the wheelchair, even when starting to drive on a thick carpet.

Reliable

LiNX follows up on the wear of both motors individually and constantly compensates. Therefore your driving comfort is guaranteed.

Intuitive interface

The LiNX interface is available on iOS and PC, allowing you to configure a chair or analyze its data quickly and easily.



Range



Separate screen

Can be positioned at different places on the chair.



Extra controller

Allows the accompanying person to take control over the wheelchair.

LiNX offers comprehensive data and statistics to monitor the performance of your wheelchair. They can be viewed in real time with a Bluetooth connection.

Tailored to your needs:

1 platform adjustable for any user

- Its modules structure and broad range of components allows you to find that perfect configuration for any user.
- All components are exchangeable and therefore the wheelchair can be adjusted over time
- Thanks to the intuitive user interface adjusting the programming is very easy.



The dongle key can be used to adjust any chair equipped with a LiNX control.



10-way switch

Fingertip control over various adjustments of the wheelchair.



Special controllers

A number of special adaptive controllers are available to find the most ideal solution for each user.

Options

Vermeiren offers a wide choice of options and controllers which are adaptable to your electronic wheelchair.

Everything is designed to improve comfort and facilitate the usage of the patient and companion. Vermeiren's solutions allow you to customize your wheelchair so that everyone can find the solution that will push the limits of mobility.



Controllers



SE64
remote controller



SE82
button controllers



SE55 - reinforced
controller



SE02 - table
controller



SE83 - mini
joystick



SE09 - Second
controller



SE04 - chin
controller



SE53 - Swing away
controller to the inside
or outside and
adjustable in height



SE53Z - Swing away
controller to the inside
or outside and
adjustable in height
with table access



SE53 A - Swing away
controller adaptive to
possible hand
deformations

Options



L04 - side supports



L58 - adjustable
headrest



SE35 - pavement
climber



B15 - plexi table



E SE18 - wheel locking system



AA2 - electrical
backrest



SE39 - electrical
seat and backrest



SE63 - softball j
oystick



SE42 - lift



BZ7 / BZ8
adjustable legrests
7: long - 8: short




BZ7e / BZ8e
electrical legrests
7: long - 8: short

Need information on spare parts or accessories ?

visit our website : www.vermeiren.com


Order forms




Videos
(product presentations,
assembly instructions...)



Exploded views
and spare parts



Measurements




390
420

225

Certificates
(CE, CERAH)
and manuals



Contacts



Compatibility table

	708 kids	Eclips X4 kids	Eclips X4 90°	708D	D200	D200 30°	Eclips	Eclips 30°	V300	V300 30°	V300 XR	V500	V500 30°	Eclips X4	Eclips AD	Sagitta	Sagita kids
B02	•	-	•	•	-	-	•	-	•	•	•	-	-	•	•	-	-
B03	•	•	•	•	•	•	•	•	•	•	•	-	-	•	•	-	-
B05	-	•	•	•	•	•	•	•	•	•	•	-	-	•	•	-	-
B06	•	•	•	•	•	•	•	•	•	•	•	-	•	•	•	•	-
BA6	•	•	•	•	-	-	•	•	-	-	•	-	-	•	•	•	-
BK6	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	•	•
BL6	•	-	-	•	-	-	•	•	-	-	•	-	-	•	-	-	-
BL7	•	•	•	•	-	-	•	•	•	•	•	•	•	•	•	-	-
BK7	•	-	•	•	-	-	•	•	•	•	•	•	•	•	•	-	-
BN6	•	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-
B07	•	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-
BZ7	•	-	•	•	-	-	•	•	-	-	-	-	-	•	•	-	-
B08	•	•	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-
BZ8	•	-	•	•	-	-	•	•	•	•	•	•	•	•	•	-	-
B09	•	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-
BZ9	•	-	•	•	-	•	•	•	-	•	-	-	-	•	•	-	-
B10	•	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-
BZ10	•	•	•	•	-	•	•	•	-	•	-	-	-	•	•	-	-
B12	•	•	•	•	•	•	•	•	•	•	•	-	-	•	•	-	-
B13	•	-	•	•	•	•	•	•	•	•	•	-	-	•	•	-	-
B14	•	-	•	•	•	•	•	•	•	•	-	•	•	•	•	-	-
B15	-	-	•	•	•	•	•	-	•	•	-	-	-	•	•	-	-
B17	-	•	•	•	-	-	-	-	-	-	-	-	-	-	•	-	-
B18	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	-	-
B19	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	•	•
B20	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B21	•	-	•	•	•	•	•	•	•	•	-	•	-	•	•	-	-
B22	•	•	•	•	•	•	-	-	•	•	•	•	•	-	•	-	-
B23	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	•	•
B24	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	-	-
B25	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	-	-
B26	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	-	-
B27	•	-	•	•	•	•	•	•	•	•	•	•	-	•	•	-	-
B28	•	•	•	•	•	•	•	-	-	-	-	-	-	•	•	-	-
B30	•	-	•	•	•	•	•	•	•	•	-	•	-	•	•	-	-
B31	•	-	•	•	•	•	•	•	•	•	•	•	•	•	•	-	-
B33	-	-	-	•	•	•	•	•	•	•	•	•	•	-	•	-	-
B37	-	•	•	•	-	-	•	•	•	•	-	-	-	•	•	-	-
B40	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	-	-
B41	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	•	•
B42	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B45	•	•	-	•	•	-	•	•	•	-	-	-	-	•	•	-	-
B46	•	-	-	•	•	-	•	•	•	-	-	-	-	•	•	-	-
	708 kids	Eclips X4 kids	Eclips X4 90°	708D	D200	D200 30°	Eclips	Eclips 30°	V300	V300 30°	V300 XR	V500	V500 30°	Eclips X4	Eclips AD	Sagitta	Sagita kids

Subject to technical modifications.

Compatibility table

	708 kids	Eclips X4 kids	Eclips X4 90°	708D	D200	D200 30°	Eclips	Eclips 30°	V300	V300 30°	V300 XR	V500	V500 30°	Eclips X4	Eclips AD	Sagitta	Sagitta kids
Bi5	-	-	-	-	-	-	-	-	-	-	-	●	●	-	-	-	-
B47	-	-	●	●	●	●	●	●	●	●	-	●	●	●	●	-	-
B52	●	-	●	●	●	●	●	-	●	●	-	●	●	●	●	-	-
B57	●	●	●	●	●	●	●	●	●	●	●	●	-	●	●	●	-
B58	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
B59	●	-	●	●	●	●	●	●	●	●	-	-	-	●	●	-	-
B63	-	●	●	●	●	●	●	●	●	●	-	-	-	●	●	-	-
B64	-	●	●	●	●	●	●	●	●	●	-	-	-	●	●	-	-
B66	-	●	●	●	●	●	●	●	●	●	-	-	-	●	●	-	-
B69	-	-	-	●	-	-	●	●	-	-	-	-	-	●	-	-	-
B73	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-
B74	●	●	●	●	●	●	●	●	●	●	-	●	●	●	●	-	-
B78	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
B79	-	-	●	●	-	-	●	●	●	●	-	-	-	●	●	-	-
B80	●	●	●	●	●	●	●	●	●	●	●	●	-	●	●	●	●
B82	-	●	-	-	-	-	-	-	-	-	●	-	-	●	-	●	●
B83	-	●	-	●	-	-	●	-	●	-	-	●	●	●	-	-	-
B84	●	●	●	●	-	-	●	●	●	●	-	●	-	●	●	-	-
B85	●	●	●	●	●	●	●	●	●	●	●	-	-	●	●	●	●
B90	-	●	●	●	●	●	●	●	●	●	-	●	●	●	●	-	-
B91	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-
B92	-	-	-	●	-	-	●	-	●	-	●	-	-	●	●	-	-
B94	-	●	●	●	●	●	●	●	●	●	●	-	-	●	●	-	-
B99	-	-	●	●	●	●	●	●	●	●	●	-	-	●	●	-	-
B101	●	-	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-
B102	-	-	●	●	-	●	-	-	-	●	-	-	-	-	-	-	-
B112	-	-	●	●	-	-	-	-	-	-	-	-	-	-	●	-	-
B115	-	-	-	●	-	-	●	-	-	-	-	-	-	●	-	-	-
AMP	-	-	●	●	-	-	●	-	●	●	●	-	-	●	●	-	-
S1	●	●	●	●	●	●	●	●	●	●	●	●	-	●	●	-	-
HEM2	-	-	●	●	●	●	●	●	●	●	-	●	-	●	●	-	-
T30	-	●	●	●	●	●	●	●	●	●	-	●	●	●	●	-	-
L04	-	●	●	●	-	-	●	●	●	●	-	●	●	●	●	-	-
L07	-	●	●	●	-	-	●	●	●	●	-	-	-	●	●	-	-
L14	-	●	●	●	-	-	●	●	●	●	-	●	●	●	●	-	-
L15	-	●	●	●	-	-	●	●	●	●	-	●	●	●	●	-	-
L24	-	●	●	●	-	-	●	●	-	-	-	-	-	●	●	-	-
L25	-	●	●	●	-	-	●	●	●	●	-	●	-	●	●	-	-
L55	-	●	●	●	-	-	●	●	●	●	-	●	●	●	●	-	-
L34	-	-	-	●	-	-	-	-	●	●	-	-	●	-	-	-	-
L35	-	-	-	●	-	-	-	-	●	●	-	-	●	-	-	-	-

● possible - not compatible

Subject to technical modifications.

Compatibility table

	Forest kids	Springer	Navix	Navix S.U.	Express	Forest 3	Forest 3 S.U.	Forest 3 Initial	Forest 3 plus	Timix	Inovys II & M.C.	Inovys II E	Serenys	V300 XXL	28	Eclips XXL
B02	-	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•
B03	-	-	-	-	•	-	-	-	-	-	-	-	•	•	•	-
B05	-	-	-	-	-	-	-	-	-	-	•	•	•	•	•	-
B06	-	•	•	-	•	•	-	•	•	-	-	-	•	•	•	-
BA6	-	•	-	-	•	•	-	-	•	-	-	-	-	-	•	-
BK6	•	•	•	-	•	•	-	-	•	•	-	-	•	-	•	•
BL6	-	-	-	-	•	-	-	-	-	-	-	-	-	-	•	-
BL7	-	-	-	-	-	-	-	-	-	-	•	•	-	-	•	-
BK7	-	-	-	-	-	-	-	-	-	-	•	•	-	-	•	-
BN6	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-
B07	-	-	-	-	•	-	-	-	-	-	-	-	-	-	•	-
BZ7	-	-	•	-	•	•	-	•	•	•	-	-	•	-	•	-
B08	•	-	-	-	•	-	-	-	-	-	-	-	-	-	•	-
BZ8	•	•	•	-	•	•	-	•	•	•	•	•	-	•	•	-
B09	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-
BZ9	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-
B10	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-
BZ10	-	-	-	•	-	-	•	-	-	•	-	-	-	•	•	-
B12	•	-	-	-	•	-	-	-	•	-	•	•	•	•	•	-
B13	•	•	-	-	•	•	-	•	-	-	-	-	•	•	•	-
B14	-	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•
B15	•	•	•	•	-	•	-	•	-	•	•	•	•	•	•	-
B17	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-
B18	•	•	•	-	•	•	•	•	-	•	•	•	•	•	•	•
B19	•	•	•	-	•	•	•	-	-	•	•	•	•	•	•	•
B20	•	•	•	-	•	•	•	•	-	•	•	•	•	•	•	•
B21	•	•	•	-	•	•	•	-	-	•	•	•	•	-	•	•
B22	•	•	-	-	•	•	-	-	-	-	•	•	•	•	•	•
B23	•	•	•	•	•	•	-	-	-	•	•	•	•	•	•	-
B24	-	-	•	-	•	•	-	-	•	•	•	•	•	-	•	-
B25	-	-	-	-	•	•	-	-	-	-	•	•	•	-	•	-
B26	-	-	-	-	•	•	-	-	-	-	•	•	•	-	•	-
B27	-	-	-	-	•	•	-	-	-	-	•	•	•	-	•	-
B28	•	•	-	-	-	•	-	-	-	-	-	-	-	-	•	•
B30	-	-	-	-	-	-	-	-	-	-	•	•	•	•	•	-
B31	•	•	•	-	•	•	-	-	•	•	-	-	•	•	•	•
B33	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B37	-	-	-	-	-	-	-	-	-	-	-	-	•	-	•	-
B40	-	-	-	-	-	-	-	-	-	-	•	•	•	•	•	•
B41	•	•	-	-	•	-	-	-	-	-	•	•	•	•	•	•
B42	•	•	•	-	•	•	-	-	-	•	•	•	•	•	•	•
B45	-	-	-	-	•	-	-	-	-	-	-	-	-	-	•	-
B46	-	-	-	-	•	-	-	-	-	-	-	-	-	-	•	-

Subject to technical modifications.

Compatibility table

	Forest kids	Springer	Navix	Navix S.U.	Express	Forest 3	Forest 3 S.U.	Forest 3 Initial	Forest 3 plus	Timix	Inovys II & M.C.	Inovys II E	Serenys	V300 XXL	28	Eclips XXL
B47	-	-	-	-	●	-	-	-	-	-	-	-	-	●	●	●
B52	-	-	●	-	-	-	-	-	-	●	●	●	●	●	●	●
B57	●	●	-	-	●	●	-	-	-	-	●	●	●	-	●	●
B58	●	●	●	●	●	●	●	-	-	●	●	●	●	●	●	●
B59	●	●	-	-	●	●	-	-	-	-	●	●	●	-	●	-
B63	-	-	-	-	●	-	-	-	-	-	●	●	●	-	●	-
B64	-	-	-	-	●	-	-	-	-	-	●	●	●	-	●	-
B66	-	-	-	-	●	-	-	-	-	-	●	●	-	●	●	●
B69	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B73	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B74	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	●
B78	●	●	●	-	●	●	●	●	-	●	●	●	●	●	●	●
B79	●	-	-	-	●	-	-	-	-	-	-	-	-	●	●	-
B80	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	●
B82	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B83	-	-	-	-	-	-	-	-	-	-	●	●	-	●	●	-
B84	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	-
B85	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	-
B90	-	-	-	-	●	-	-	-	-	-	●	●	-	●	●	●
B91	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	-
B92	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-
B94	-	-	-	-	●	-	-	-	-	-	●	●	●	●	●	-
B99	-	-	-	-	●	●	-	-	-	-	-	-	-	-	●	-
B101	-	●	-	-	●	●	-	-	-	-	-	-	●	●	●	-
B102	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-
B112	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-
B115	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	-
AMP	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	-
S1	-	-	-	-	-	-	-	-	-	-	●	●	-	●	●	-
HEM2	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	-
T30	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	-
L04	●	-	-	-	-	●	-	-	●	-	●	●	●	●	●	●
L07	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	-
L14	●	-	-	-	●	-	-	-	-	-	-	-	-	●	●	●
L15	●	-	-	-	●	-	-	-	-	-	-	-	-	●	●	●
L24	-	-	-	-	●	-	-	-	-	-	-	-	-	-	●	-
L25	-	-	-	-	●	-	-	-	-	-	-	-	-	-	●	-
L55	●	●	-	●	●	●	●	●	●	-	-	-	●	●	●	●
L34	-	-	●	●	●	-	●	-	●	●	●	●	●	●	●	-
L35	-	-	●	●	●	-	●	-	●	●	●	●	●	●	●	-

● possible - not compatible

Subject to technical modifications.



electro

Scooters



Antares 3 - Antares 4.....	68
Venus 3V - 4V - Venus 4 Sport	69
Sedna - Sedna Premium	70
Luggie - Luggie Super	71
Ceres 3 - 4 Deluxe	72
Ceres S.E.....	73
Eris	74
Carpo 3 - Carpo 4 Limited	75
Carpo 2 Standard - Eco - Sport	76

Carpo 2 XD S.E.	78
Carpo 2 S.E.....	79
Mercurius 4 Ltd.....	80
Options	82
ScooterPac	83
Ramp kits	84
Telescopic / Custom ramps	85

Antares 3

FEATURES
Disassemble in 1,2,3,4...5 parts



blue



6
km/h

The scooters Antares 3 and Antares 4 are perfect for indoor usage.

Easy to store and to move around and can be disassembled in smaller parts for easy transport.

Adjustable steering

Basket

Stability wheels

Swivelchair adjustable in height

Full tires

Battery with firm handle to easily remove the battery

Antares 4



blue



6
km/h



Easy to disassemble for transport

Adjustable steering

Basket

Swivelchair adjustable in height

Full tires


Battery with firm handle to easily remove the battery



km/h



km



A3

6

A4

6

12

12

8°/14%

8°/14%


54

54

1950

2790

Options and accessories pag. 82



A3

480

A4

480

980

980

965

970

200 x 50 / 3

200 x 50 / 4

310-360

310-360

41

43

115

115

12V / 12 Ah

12V / 12 Ah

270 W

270 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Venus 3V & 4V

FEATURES

Disassemble in 1,2,3,4...5 parts



blue



6
km/h



optional
20A batteries
available!

Adjustable steering

Basket



Swivel chair
adjustable in height

optional
20A batteries

Full tires

Venus 4 sport



champagne



6
km/h

A compact 4-wheel scooter for indoor and outdoor usage. Comfortable and practical in day to day usage and can be disassembled in 5 parts without extra tools.



km/h



km



Frontlight


Basket



Full tires

Swivel chair
adjustable
in height

Options and accessories pag. 82



V3V	510	900	1090	200 x 50 / 3	525-575	43	130	2 x 12V / 12 Ah	270 W
V4V	510	900	1070	200 x 50 / 4	525-575	46	130	2 x 12V / 12 Ah	270 W
V4S	550	900	1040	220 x 80 PU/4	500-550	55	115	2 x 12V / 20 Ah	270 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

VENUS 3V - 4V

Sedna (Premium)

FOLDABLE SCOOTERS



blue
(standard)



red
(premium)



Sedna scooters are perfect to carry along as they are foldable in a few minutes. They are designed with a removable lithium battery with anti-theft security and inflatable tires. The Premium version (red) is equipped with a swivel seat and armrests. The standard version (blue) only offers them as an option.



STEP 1



STEP 2



Option: Transport net



STEP 3

FEATURES

Easy foldable scooter

Folable backrest

Adjustable column

Foldable and pivoting chair
(Premium version only)

Inflatable tires

Anti-slip mat

Lithium battery
with anti-theft lock

STEP 4



8



20



6°/10%




30



820



530



840



870*



2.80/2.50 - 4 air(2)
200 x 50 air(2)



540



27,3

*with armrests
and swivel chair



115



24V / 11,5Ah



270W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Luggie (Super)

FOLDABLE SCOOTERS



blue
(Super)



champagne
(Standard)



6
km/h

The Luggie is the ideal transportable scooter. Convenient for short distances, foldable in minutes and transportable as a carry-on suitcase. Very convenient to drive in shops or restaurants.



Easily foldable for car transport




Extendable handles
(on the STANDARD version)




Armrests (optional on the
STANDARD version and
standard on the SUPER version)



LUGGIE - LUGGIE SUPER



Standard	6	20	6°	25	1038
Super	6	18	6°	25	1050



Standard	455	860	982	6" front(2) 7" rear(2) PU	120	25,4	114	8,5Ah
Super	525	860	989		120	30,5	164	10,5Ah

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Ceres 3 - 4 Deluxe



champagne

model CERES 3 - 4 deluxe

Carpo's little brothers are comfortable scooters for your daily usage.

FEATURES

Comfort and style

Electronically programmable



Adjustable swivel seat
with removable headrest
and adjustable armrests

- Mirror
- Left and right
indicators

Shopping basket

Inflatable tires



Delta steering

Bumper

Inflatable tires

Adjustable seat

Anti-slip floor mat



km/h



km



6°/10%
6°/10%



100
100




3020
2940

C3
C4

12
12

35
35

Options and accessories pag. 82



C3
C4

610
610

1160
1160

1270
1310

330 x 100
330 x 100

445-520
445-520

92
94,7

135
135

2 x 12V / 50 Ah
2 x 12V / 50 Ah

470 W
470 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Ceres S.E.

FEATURES

The most compact of outdoor scooters



silver grey



red

The Ceres S.E. (Special Edition) combines design and quality. A simple manual controller, new LED lighting and sports bumper for safety.



Mirrors (2)



Manual control panel with big and visible icons

Great comfort with adjustable seat, headrest and armrests

Transport basket

Front and rear lighting


Bright indicators

Angle adjustable steering column

Anti-slip floor mat

Lateral reflectors for better visibility

Bumper



10

35

6°/10%

115

2800

15


48

30

Options et accessoires page 82



620



1280



1300



4.00-5 air (4)



175



106



140



2 x 12V / 38Ah
2 x 12V / 50Ah



450W


Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

CERES S.E.

Eris

INDOOR / OUTDOOR SCOOTER

The Eris model is an indoor / outdoor scooter ideal for small distances between shops. Equipped with indicators and LED lighting. Includes seat adjustment and shopping cart.



Removable and in width adjustable armrests

Height adjustable swivel seat

LED lights

Adjustable steering column

Anti-slip floor mat


Lateral reflectors for a better visibility



6
10



35-45
30-37



6°



50




2900

Options et accessoires page 82



570



1150



1185



3.00-4
air (4)



145



75



140



2 x 12V / 38Ah
2 x 12V / 50Ah



320W
320W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Carpo 3 & 4

FEATURES

Speed and lasting pleasure



brown

Electronically programmable



Model Carpo 3 Limited edition

Powerful outdoor comfortable 4-wheel scooter. Electronic programmable and a reliable partner for outdoor trips.

optional
80A batteries
available!

Adjustable swivel seat
with removable headrest
and adjustable armrests

Lighting front and
rear and indicators

Bumper

Inflatable tires

Extra strong
shopping basket

Model Carpo 4 Limited edition

optional
80A batteries
available!

Powerful and comfortable 4 wheel outdoor scooter. The Carpo 4 Limited edition is electronically programmable and standard equipped with a bright LCD with various indications.

Lighting front and rear
and indicators

Bumper

Inflatable tires

Adjustable seat



Also available as
Model Carpo 4
in champagne color



champagne



km/h



km



9°/15%
9°/15%




100
100



3030
3160

available
XD version
For hilly regions

Options and accessories pag. 82



C3 660

1200

1400

4.00-4/3.00-8

460/480/505/530

127

150

2 x 12V / 75 Ah

750 W

C4 660

1350

1450

3.00-8

455/480/505/530

131

150

2 x 12V / 75 Ah

750 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

CARPO 3 - 4

Carpo 2 Standard



champagne

This strong and comfortable 4-wheel scooter is ideal for outdoor usage. It can be programmed electronically and provides you excellent driving pleasure.

Electronically programmable

Optional



brown



black



blue



red



Electronic display



Extra strong shopping basket

Comfortable adjustable seat and headrest

Battery 75A (optional 90A)

Front and rear lighting and indicators

Bumper

Anti-slip floor mat

Inflatable tires



10/15/18



34-65



6°/10%



100




2000

Options and accessories pag. 82



660



1420



1648



100x65-9
air- (4)



490-565



144,6



150



2 x 12V / 90 Ah
2 x 12V / 75 Ah
2 x 12V / 90 Ah



450 W
600 W
750 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Carpo 2 ECO

Electronically programmable



blue

Optional



brown



black



champagne



red



Analog display



Optional: Spare tyre



Optional: Luggage box

Carpo 2 Sport

Electronically programmable



brown

Optional



blue



black



champagne



red



Electronic display



Phone compartment
USB connector

Carpo 2 XD S.E.

FEATURES

Special for hilly regions

OUTDOOR SCOOTER



silver grey



red



Programmable electronics



- USB connector
- Phone compartment

Extra comfort with comfortable seat, headrest and armrests

Solid, portable shopping basket

Armrests adjustable in height and width

Batteries 90A

Front and rear light and indicators

Anti-slip floor mat

Inflatable tires & black rims

Bumper



13



57



6%/10%



100




4020

Options and accessories pag. 82



660



1390



1535



3.00-8
(4)



210



144



150



2 x 12V / 90 Ah



750 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Carpo 2 S.E.

FEATURES

Long distance drivers

OUTDOOR SCOOTER



silver
grey



Solid,
portable
shopping basket



Programmable electronics

Extra comfort with
comfortable seat,
headrest and
armrests

USB connector

Phone
compartment

Armrests
adjustable
in height
and width

Front and rear
light and
indicators

Batteries 90A

Anti-slip
floor mat

Bumper

Inflatable tires
& black rims



10/15



60



6°/10%



100




4020

Options et accessoires page 70-71



660



1390



1535



3.00-8
air (4)



210



144



150



2 x 12V / 90 Ah



750 W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

CARPO 2 S.E.

Mercurius 4 LTD.

FEATURES
For long distance drivers

OUTDOOR SCOOTER



champagne



Programmable electronics

Extra comfort with comfortable seat, headrest and armrests

Mirrors

Solid, portable shopping basket

LED lights front and back and indicators

Double USB connector

Armrests adjustable in height and width

Batteries 75A (option 80A)

Suspension

Inflatable tires & sports rims

Lateral reflectors for a better visibility



15



40



8°/14%



100




3740

Options and accessories pag. 82



680



1355



1430



3.00-8 (4)



165



116,75



150



2 x 12V / 75 Ah



1000W

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Options

Accessories : Vermeiren offers a wide selection of accessories for your scooter.



Crutch holder (SW31 - SH31)



Canister holder (B83)



ScooterPac - rain protection



Protective cover (SH08)



Mirror (x1)



Rollator holder (SH06 - SW06)



Luggage box (SW03)

options	1. Crutch holder	2. Canister holder	3. Scooterpac	4. Protective cover	5. Mirror	6. Rollator holder	7. Luggage box
1. Crutch holder				OK	OK		
2. Canister holder				OK	OK	OK	OK
3. Scooterpac					OK *		
4. Protective cover	OK	OK	OK *		OK	OK *	OK
5. Mirror	OK	OK		OK		OK	OK
6. Rollator holder		OK		OK	OK		
7. Luggage box		OK		OK	OK		

* The mirror can be mounted but not used while driving when the ScooterPac is unfolded.

options	Anlares 4	Venus 3V/4V	Venus 4 sport	Ceres 3 Ceres 4 S.E.	Carpo 4	Carpo 2 ECO	Carpo 2	Carpo 2 S.E.	Mercurius 4 LTD
1. SW(SH)B31	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2. B83	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3. Scooterpac				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
4. SH08				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5. SH(SW)74	<input type="checkbox"/>		<input type="checkbox"/>	S	S	S	S	S	S
6. SH(SW)06			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7. SV07						<input type="checkbox"/>	<input type="checkbox"/>		
8. SW03			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

S standard

☐ compatible option

The world's first fold-away mobility scooter canopy



Folds out in seconds
for protection from
any weather conditions

**Enjoy independence with this
easy-to-use canopy and go out
on a whim, without any help.**

Looks like rain?

**Fold it out from behind the seat
to form a top and side canopy.**

Turned out sunny?

**Fold the canopy back in seconds
for a breezy ride.**



SCOOTERPAC

Some of ScooterPac's unique features



Speedy

No long waits for a canopy to be made, and with the added benefit of it fitting virtually every scooter - instantly!



Versatile

From the largest 9mph scooters right through to the smallest car-transportable, the ScooterPac's clever design means one size fits all!



Innovative

The ground breaking patent-pending technology provides protection from the rain and lets you soak up the sun (when it's out!)



Robust

Aluminium construction for long term durability and no chance of rust! And-state-of-the-art cover material too!

Quickramps



These ramps are well suited for small obstacles such as interior thresholds for walking with a rollator


Quickramp is an ideal ramp for obstacles of 10 to 16 mm.

Available in 4 different heights: 10 mm, 12 mm, 14 mm and 16 mm and made of the same solid polyethylene material as the larger ramps.


Ramps

BUILDING KITS

model n°1




model n°3



model n°2






model n°4



Which ramp?

Determine the height to cross and then select the proper ramp kit. From 10 mm to 150 mm, this modular system allows you to fine tune the kit to your desired height and support up to **200 kg**.



Model

Kit 1
Kit 2
Kit 3
Kit 4

260
510
760
1010

750
750
750
750

10 - 40
45 - 75
75 - 115
115 - 150

1,3
4,3
9,2
15,8

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Ramps



2 m
model



Telescopic aluminum ramps with non-slip surface for maximum safety.






3 m
model



Customized Vermeiren products

Vermeiren allows you to create your own customized ramp system.
Kindly contact our sales department.



Model				
691	1060 - 1900	250	50	7 x 2
691	1660 - 3130	250	50	10 x 2

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Comfort

Comfort chairs

Inovys II - Inovys II-E	88
Serenys	90
EZ1	92

Shell chairs

Charme - Charme E	94
Regina	95
Altitude	96
Coraille	97
378	98

Inovys II

MULTI COMFORT CHAIR

- Adjustment of the armrests in length and height
- Adjustment of the seating depth
- Adjustment of the footrests in width
- Possible to lower the chair without changing its original wheels



C55
inox

Headrest adjustable in angle,
height and depth with
rotation and translation

Armrests can be removed.
Adjustable in height, length
and width

Memory foam
seat cushion

Pushing handles
adjustable in height



Complete tilting for total relaxation



Easy to disassemble for transport or storage




*Optional L04
Side supports*



Optional B12 - Table

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



630	390-440	220-320	460-495-530	1060	550-585	0°-21°	0°-54°	460-520	1010-1050	135	36
680	440-490										
730	500-540										

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Inovys II-E

FEATURES

Electrical adjustments
for a greater comfort

MULTI COMFORT CHAIR WITH ELECTRICAL SEAT AND BACKREST

The Inovys II-E is a highly comfortable chair with multiple adjustment possibilities and an independent battery.



C55
inox



C60
brown

Headrest adjustable
in angle, height
and depth



Removable armrests.
Adjustable in height
and width



Electrical adjustment of the seat
and the backrest


Push handles
adjustable
in height



Memory foam
seat cushion

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



630
680
730

390-440
440-490
500-540

220-320

460/495/530

1060

550-585

0°-21°

0°-54°

460-520

1010-1050

135

43,9

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

INOVS II-E

Very comfortable chair for prolonged seating.
Equipped with memory foam cushions, wraparound armrests, a deformable headrest and legrests with retractable cushions facilitating transfers.
Smooth tilting back and forth due to two springs.
Possibility to lower the seating to 450 mm (with standard cushion) without changing the wheel size.



C20
blue



C29
grey



C41
green

FEATURES

Foldable and easy to transport

Headrest adjustable in angle, height and depth



Height ground /seat adjustable
in 3 positions

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37)
and check the compatibility table on page 62 - 65.



630
690
750

390
450
510

220-290

420-520
(front)
380-530
(rear)

1020

500/600

19°/113°

-5°/50°

-2°/15°

450-500

1170


120



*Comfortable, height adjustable armrests.
Retractable for easy transfers.*




*Long backrest: 600 mm (order only)
Short backrest: 500 mm (standard version)*



*Double gas cylinders providing assistance
when the patient wants to get up.*

The legrests are fully adjustable and the footrests with inclination will give an ideal angle for the feet at any given time.



*This chair is the most comfortable one on the market for prolonged seating due to its memory-foam seat cushion.
Foam density: 41 kg/m³.*



Seat angle and seat height adjustable

1. The legrests are adjustable in height.
2. The calf supports are adjustable in height, width and depth.
3. The calf supports are collapsible to the back to ease transfers.



EZ1

MULTI COMFORT CHAIR



C29
grey

Headrest with several adjustments



The **EZ1** wheelchair model offers you a wide variety in options and possibilities. For the ultimate comfort in combination with lasting materials we are pleased to introduce the EZ1 to you.

Armrests removable and adjustable in height



Push handles adjustable in height

Anti-tipping wheels



The footplates have their own unique adjustment variations which enables the angle of the footplate to be adjusted to achieve perfect comfort.



L04 - Side support



Removable and adjustable armrests, seat adjustable in depth



B12 - Wooden table



Model with headrest, side supports and table



Adjustable width

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.

620	390											
650	420											
680	450											
710	480											
740	510											
770	540											
800	570											
		190-290	450-600	1020	500/600	19°-113°	30°	-5°-15°	430-520	1185	150	
					560							

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Charme

SHELL CHAIRS

The Charme shell chair is the result of design, comfort and aesthetics for people with a weakened muscle tone. The chair allows the user to easily position himself in a resting position and keep his autonomy. It integrates in the patients home due to its looks and design and adds a therapeutic value.

Coordinated tilt of seat and legrest for the perfect positioning without complicated adjustments.

Two coloured upholsteries for a perfect match with existing interiors are available.

STANDARD VERSION



kaki beige
fabric



artificial leather
taupe brown fabric

DELUXE VERSION



artificial leather
black



artificial
leather chocolat

FEATURES

Manual or electrical



Manual version



378 - Table on wheels



780



380
400
420
440
480
270
(cushion 60 mm)
250
(cushion 80 mm)
230
(cushion 100 mm)



450-475
470-495
490-515



1220-1245
1220-1245
1260-1285



750
770
790



-10/+30°



450-475
470-495
490-515



700-1280



25
29

(electrical version)



120

Regina

DESIGN SHELL CHAIR

The Regina chair was designed for people whose muscular tones have weakened. The electrical inclination, operated by a remote control, allows the user to position him or herself in a nice resting posture and still maintain full autonomy.

This particularly neat design takes the chair away from the medical environment and integrates it perfectly into your house.

black skai
black fabric

chocolat skai
brown fabric

- Power tilt by remote control
- Shaped back with lumbar support
- Concealed chassis with an elegant design



Comfortable legrest with automatic lengthening



378 - Table on wheels

Removable headrest

Plastic wheels with brakes



700
925 (avec table)

380
400
420
440
480

230-250-270

495-520
(sans coussin)

1330-1355

665

500

-3° / +28°

35

120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

REGINA

Altitude


FEATURES

Easy maintenance

FUNCTIONAL SHELL CHAIR

The Altitude, in its basic form, is equipped with an 80 mm thick seating cushion (backheight: 770 mm). The height of the back and armrests are depending on the thickness of the seating cushion (60-80-100 mm).

Optional anti-decubitus cushion available.



Optional B52
Serum holder



Optional separated footplates




378 - Table on wheels

Option: **B12** - Removable
table. Collapsible and
adjustable in depth
(for chairs with
seatwidth 380 - 480 mm)

Footrest in one piece
mounted on a spring
in order to stay upright
during a transfer

Also available: Optional anti-decubitus
Dartex cushion (80 mm)



670	380 400 420 440 460 480 500 520	230 (cushion 60mm) 250 (cushion 80mm) 270 (cushion 100mm)	450.470.490 510.520.530 540.550.560	1190 1215 1240	750 770 790	-4/36°	460	1120	36	120
790 (XXL) 790 (XXL)	600	230	520.540.560	1240	780		500	1120	37	150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

FUNCTIONAL SHELL CHAIR

Shell chair for indoor / outdoor usage
with 4 wheels: 2 big rear wheels and 2 front wheels.

Optional anti-decubitus cushion available.



dartex
green



dartex
purple



dartex
praline



dartex
light purple



Optional B52
Serum holder



378 - Table on wheels



Optional separated footplates



Also available: Optional anti-decubitus
Dartex cushion (80 mm)



670	380 400 420 440 480 550 600	230 (cushion 60mm) 250 (cushion 80mm) 270 (cushion 100mm)	450.470.490 510.520.530 540.550.560	1190 1215 1240	750 770 790	-4/36°	460	1120	36 37	120 150
790 (XXL) 790 (XXL)		230	520.540.560	1240	780		500	1120		

378

RECLINING TABLE

This wooden table with metal frame can be used interchangeably with several Vermeiren products such as a geriatric chairs, relaxing chair,... It can be set in several tilt positions and each wheel has its individual brake. The height of the table is also adjustable according to your needs.



brown



600



507



957



615-805



0°-90°



11,15



10

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

beds & Homecare

Beds

Luna 2 - Luna 2 basic	100
Luna 2 ultralow	101
Interval	102
Alzeis	103
Pratic	104

Patient lifts

Eagle	105
Albatros 1 & 2	106

Geriatric chairs

Normandie	108
Languedoc	109
Provence	109
Alesia	110

Lift chairs

Ontario 1 - 2 motors	112
Montreal	113

Luna 2 & Luna 2 basic

This standard bed can also be used as a bed at low height because its frame can be positioned at 2 different heights (1-2). The lowest height can be 330 or 440 mm and the maximum load is 175 kg (mattress & patient).

Available with metal barriers

- Electric head part
- Electric knee buckle, manual lower leg part
- 2 pairs of side barriers in wood positionable in two different heights
- Self-raiser to be positioned left or right
- 4 directional wheels with brakes




beige



Luna 2 basic




Version with metal side barriers




The frame can be positioned at two different heights to lower the bed.




Highest position



Transport kit



Luna 2



900x2000x120
max. 20 kg

1025

2170

335 / 400

745 / 810

31,7°

83°

-14,8°

140

98,15

1015
(Luna 2 basic)

83,70
(Luna 2 basic)

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Luna 2 Ultra Low

EXTRA LOW POSITION



beige

The Luna UL 2 is a 4 section electric profiling bed that is ideal for those confused or disorientated patients who may be prone to falling out of bed. With the variable height position in its lowest setting and with a soft foam roll-out floor pad in place the dangers of injury from unscheduled bed exit can be minimised.


- The minimum height from the ground is 205 mm
- Electrical adjustable with the included controller
- Mattress size 900 mm x 2000 mm
- Controller with lock function
- Electrical head adjustment
- Electrical knee buckle, the foot part can be adjusted manually
- Two pair of wooden siderails
- Two pair of wooden siderails
- Stylish front and rear wooden boards
- Easy to (dis)assemble



Anti-trendelenburg



Highest position



900x2000x120
max. 20 kg

1025

2180

205

615

si

31,7°

83°

-14,8°

140

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Interval

CLASSIC BED

The Interval bed has an electrically adjustable head support. By simply pressing a button, the head support moves into a more vertical position, enabling you to read, eat or watch television. The footrest can be placed in the correct position manually.

This high/low bed has 4 wheels, each with a separate brake, therefore guaranteeing stability. The strong and flexible coated baseboard is easy to maintain. The bed can be assembled and disassembled quickly and without tools.



Optional centralized brake available




beige

Options

- Bed with electrical knee buckle
- Table
- Bedside table
- Cupboard - table - chair
- Cover to transport the wooden panels
- Wooden base frame
- Centralized brake



kg



900x2000

1030

2080

370

805

no

31,7°

70°

30°

150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Alzeis

ROLL-OUT BED

Security, practicality and esthetics. Low height bed for patients with Alzheimer's disease. Very low minimum height and maximum height up to 720 mm.

Compatible with most patient lifts on the market.



white

Options

- 3 functions electrical bed with knee buckle
- Protective transport cover for the wooden panels
- Serum holder
- Louis Philippe wooden panels
- Remote control support
- Independent serum holder



Optional cradle for the remote




Available with metal side barriers



ALZEIS



kg



900x2000

1030

2050

200

720

no

30°

65°

16°


140

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Pratic

PRACTICAL COMFORT

The Pratic bed is equipped with a head-rest with electronic controls. A simple touch on a keypad enables you to raise it for reading, watching television or to take your meals. The legrest can also be adjusted to any desired position.



Delivered with a transport kit



Wooden side barriers



white

Adjustable self-raising device with serum holder. Possible to install on the left or right side of the bed.

Options


- Centralized brake
- Louis-Philippe panels
- Independent serum holder



*Optional Louis-Philippe panels
light oak or dark oak.*

Exists with 2 or 3 functions
Wooden or metal barriers

Next to the side rails is an additional rail to keep the mattress in its place.



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eagle 620/625

PATIENT LIFT

The Vermeiren patient lift was designed to bring a real comfort to the patient during the transfers but also for the people who are taking care of the patient.

The patient lift Eagle 620 and Eagle 625 are equipped with the Linak HB7X engine feeding from a battery with an

autonomy of approximately 45 lifts.



Users weight
Eagle 620 : 150 kg
Eagle 625 : 175 kg


Foldable model: easy storage

Slings

SANITARY

Available in
M - L - XL (sizes in mm)


	A	B	C	D
M	740	220	580	410
L	780	260	660	470
XL	860	310	710	500



COMFORT

Available in
M - L - XL (sizes in mm)


	A	B	C	D	E
M	750	350	470	270	400
L	820	380	490	400	400
XL	900	430	590	400	400



STANDARD

Available in
M - L - XL (sizes in mm)

	A	B	C	D
M	740	370	430	400
L	790	420	480	430
XL	850	480	520	450



- ✓ **Design**
- ✓ **Functionality**
- ✓ **Safety**
- ✓ **Reliability**

Quality wheels of 100 mm ensure that this patient lift can be used on almost any underground.



Model

620
625

1740
1840

1115
1290

550
680

1035
1225

728-1740
766-1840

480
480

37
42

150
175

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Albatros 1

PATIENT RAISER

Compact patient hoist with easy usage.
The Albatros model is foldable in a few seconds for easy storage and transport.




C29
grey

Slings

Sling sizes

Available in
M - XL

M	XL
40 - 120 kg	100 - 200 kg



Integrated charger



High quality wheels
with separated brakes



1046



584



945



1040 - 1739



500



44,85



150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Albatros 2

PATIENT RAISER

Compact patient hoist with easy usage. The Albatros 2 model is foldable in a few seconds for easy storage and transport. The Albatros 2 helps to stabilize a person in a verticalised position.



C29
grey



Maintain a vertical position with ease



High quality wheels
with separated brakes



1115



584



945



1380 - 1687



500



50,6



150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

ALBATROS 2

Normandie

GERIATRIC CHAIR

The Normandie is a geriatric chair with a reclinable backrest up to 45°. Equipped with a legrest that reclines automatically with the inclination of the backrest.



C32
yellow



C33
light
purple



C34
orange



C35
light
green



C36
dark
red



C37
night
blue



C38
brown



C39
deep
green



C40
blue
fabric



C41
red
fabric

OPTIONS



Headrest



Extra wide armrest
movable to the side



Adjustable in height
(from 508 to 568 mm)




Version with wheels



B12 - Table



378 - Table on wheels



Also
available
as electrical
version.



660



540



0-170



510



1200




750



440



13°-45°



780



24,85



130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Languedoc

GERIATRIC CHAIRS

PROVENCE : Geriatric chair with a fixed back and height adjustable armrests.

LANGUEDOC : Geriatric chair with a reclining back (45°) and height adjustable armrests.



C32
yellow



C33
light purple



C34
orange



C35
light green



C36
dark red



C37
night blue



C38
brown



C39
deep green



C40
blue fabric



C41
red fabric

RECLINABLE BACK 45°



OPTIONS



378 - Table on wheels



Retractable footrest



Reclining footrest
500 x 500 x 450 mm



B12 - Table



Headrest



Serum holder

Provence

FIXED BACK

Provence B17

B17 - TOILET SEAT

ALL MODELS

are originally equipped with height adjustable armrests



Languedoc
Provence

660

540

0-170

510

1200

750

440

13°-45°

780

23

130

660

540

0-170

510

1150

750

440

10°/21°/32°

780

22

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

PROVENCE

Alesia

GERIATRIC CHAIR

This geriatric chair is equipped with 125 mm (Ø) wheels for easing mobility. The backrest has been kept low to make it easy to wash the head of the patient. The backrest is reinforced in the lower back and is tiltable by the user or third party because the controls are connected. The seat can be tilted forward and backwards and the footrest palette is spring-mounted and retractable.

The armrests are adjustable in height and the chair is available with or without a legrest. The brake is centralized and the rear wheels are directional.

There are three operating pedals: the blue for tilting the chair, red for the brakes and green to block the directional wheels.



C32
yellow



C33
light
purple



C34
orange



C35
light
green



C36
dark
red



C37
night
blue



C38
brown



C39
deep
green



C40
blue
fabric



C41
red
fabric

OPTIONS

- Serum holder
- Adjustable and removable headrest
 - Table
- document holder
- Canister holder
- Urine bag holder



700



540




0-270



530



1330



800



12°



30°




480



800*
1330**



33



140

(*backrest and legrest not inclined. **backrest and legrest inclined.)

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

More comfortable seat




New adjustable footplate
(3 different height positions)



New optional wheels,
more practical and redesigned



Shortened armrest tube to
prevent contact with
the floor when the armrest
is at its minimum height



Newly designed pedal
brake for a more pleasant use



Handle repositioned for easy
access, especially when
the seat is tilted



Ontario I

LIFT CHAIR

The Ontario lift chair offers you optimal comfort in your own house. With a simple push on a button you can relax, watch tv or have a little nap in your Ontario lift chair.

With the Ontario 2 motors, you can separately control the backrest and the legrest.




Ontario II



Optional armrest covers



378 - Table on wheels



Ontario	800	465	150	500	1160	840	450	695	49,64	150
Ontario II	845	480	140	500	1145	960	530	670	53,34	150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Montreal

wall model

LIFT CHAIR

The Montreal can be placed very close to the wall (100 mm) and still recline completely.



skai
noir



skai
taupe



tissu
beige



378 - Table on wheels

The Montreal model will offer you the same possibilities as the Ontario. Make your choice out of different colors and textures for an optimal match with your existing interior. The footrest inclines first making it easy to watch television without having to lay back completely.



MONTREAL




850



485




150




500




1140



830



510



700



53,98



150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



x-tra

Bariatric products

Manual wheelchairs

V-series XXL	116
28	117
Eclips XXL	118
Forest 3 Plus	119

Geriatric chairs

Normandie - Languedoc XXL	120
---------------------------------	-----

Shell chairs

Altitude XXL	121
Coraille XXL	121

Commode chairs

9062 XXL	122
----------------	-----

Beds

Interval XXL 1200 - 1400 mm	122
-----------------------------------	-----

V-series XXL

XXL WHEELCHAIR

These extra strong wheelchairs with double cross can support a maximum weight of 170 kg.

V100 XXL - standard - steel frame

V200 XXL, V300 XXL - lightweight - aluminium frame



C59
champagne



C80
blue

170 KG

Metal pushbar

Model V300 XXL
with adjustable
backrest

Model V200 XXL

Double cross for
additional strength



Model V100 XXL with a fixed back and seat




Model V200 XXL



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



	740	770	810	375	530	560	600	220-240	470-530	380	890-950	490	1060	21,4	170
V100XXL	740	770	810	375	530	560	600	220-240	470-530	380	890-950	490	1060	21,4	170
V200XXL	740	770	810	385	530	560	600	220-240	440-530	400	820-1000	490	1050	20,1	170
V300XXL	740	770	810	350	530	560	600	220-240	440-530	400	840-930	490-530	1060	21	170

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

XXL WHEELCHAIR

150 KG

This chair with double cross and fixed back can support a weight up to 150 kg.

Equipped with retractable and foldable footrests, desktype armrests that allow easy access to a table and comfortable long armrest padding.



C20
blue



Anatomic version (optional)
Max. users weight: 150 kg



Seat and backrest in black nylon

Double cross for additional strength

Small side wheels on the footplates will protect walls and furniture while moving around

Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



560	290	390							18	150
610	290	440							18	
670	310	500	230	500	420	920	430	1110	19	
750	350	550							20	

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips XXL

XXL WHEELCHAIR

Aluminum foldable wheelchair with double cross. The Eclips XXL provides a natural comfort with its anatomical seat and backrest. This wheelchair, equipped with reinforced rear wheels, can support a maximum weight of 200 kg.

Easy transfers thanks to the foldable aluminum footplate in one piece.



C27
blue

200KG

B74 - Optional drumbrakes

Anatomic
seat and backrest

Adjustable
heel strap

Reinforced clothing protectors

Foldable, one piece, footplate



B78 - Anti-tipping



B74 - Brakes



Wide range of accessories and options.

Refer to the accessories page (page 36 - 37) and check the compatibility table on page 62 - 65.



850 900 950 1000	375	600 650 700 750	200	530	620	1140	530/540/550/560	1140	38,5 39 40 40,5	200

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Forest 3 Plus

ELECTRONIC XXL WHEELCHAIR

The Forest 3 Plus is a wheelchair that covers most of the features of a standard Forest 3 but specially designed for larger sizes.



C71
blue



C72
red



C73
green



C76
white



C77
black



C78
blue

250KG

FEATURES

Special Forest 3 for wides sizes



All new **LiNX** controller for easier movement and access

LiNX ACCESS Wheelchair Control System

LE Controller System :

- Ideal for light weight, indoor wheelchairs
- Easy to use
- Simple
- Reliable



2x12V-80Ah 2x350W 6 km/h 49 km 9° 52 1840
10 km/h 49 km

Wide range of accessories and options. Refer to the accessories page (page 60) and check the compatibility table on page 62 - 65.



705 500-550 335-380 450 450/500/550 590 250-399 1145-1163 134,35 250
760 600-650 335-380 450 450/500/550 590 250-399 1145-1163 134,35 250

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

FOREST 3 PLUS

Normandie XXL

Languedoc XXL

MAX
users weight
350 kg

See page 109 for the options

GERIATRIC XXL CHAIR

The Normandie is a geriatric chair with a reclinable backrest up to 45°. Equipped with a legrest that reclines automatically with the inclination of the backrest.



Languedoc XXL



C32
yellow



C33
light
purple



C34
orange



C35
light
green



C36
dark
red



C37
night
blue



C38
brown



C39
deep
green



C40
blue
fabric



C41
red
fabric



Armrests
adjustable
in height

Reclinable back 45°
and legrest will recline
automatically with the
back



890

770

0-220

510

1200

760

480

13°-40°

780

37,75

350

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Altitude - Coraille XXL

GERIATRIC XXL CHAIR

Shell chair Altitude and Coraille XXL
Seat width 550 mm or 600 mm for
the overall width of 790 mm
Maximum users weight: 150 kg

MAX
users weight
150 kg



dartex
green



dartex
purple



dartex
praline



dartex
light purple



378 - Table on wheels



5 points safety
harness

Only available with
separated footrests

GERIATRIC XXL



790	550	230	520.540.560	1240	1120	500	-4°/36°	780	37	150
790	600	230	520.540.560	1240	1120	500	-4°/36°	780	37	150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

9062 XXL

XXL TOILET CHAIR

XXL commode with a seat width of 600 mm and max. users weight of 175 kg.



C11
black

max. users weight
175 kg



Bucket with handle



780 600 510 540 400 15 175

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Interval XXL



Metal side rails version



max. users weight
240 kg


The 3-functions Interval XXL bed has electrical adjustments and is available in 1200 and 1400 mm.

Delivered in a white frame with a grid mattress base. The legpart is provided with an electrical adjustable knee-buckle whereas the footpart can be adjusted manually.

Max. users weight: 240 kg.



Delivered on a transportset.
Installed and ready for usage in a few minutes.



1200x2000 1330 1530 2080 370 805 - 28,5° 70,5° -20,5° 240

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

cycling

Freedom



Freedom	124
Aqua Safari Happy Sporty	125

e-Bike	126
Vintage - Lagoon	127
Options	128

FREEDOM

BLACK WHITE

Model FREEDOM.

The Freedom is a sturdy tricycle ideal for demanding adolescents available in black and white. A perfect tricycle for strolling through countryside..

Range of options available
at page 128.

14 years
and
above



Also available as
electrical version



Inside leg
from 68 to 88 cm

- The steering wheel (different steerings available) is adjustable in height and angle.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- V brakes at the front with parking lock. Optional disc brake for the rear wheels available.
- Big metal transport tray at the rear.



total
length

1599



total
width

670



saddle
height

709-950



distance
saddle-pedal

586-827



step-in
height

174



height
steering

986-1174



wheels

20" x 1,75



total
weight

25,5



users
weight

100

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Aqua

3 years
and
above



Inside leg
from 30 to 46 cm

Model AQUA. Decorative figures and bright, shiny colours form the main attraction for young children. It has a bell with a cute shape and a rear tray, ideal to transport toys around.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- The steering wheel (different steerings available - see page 128) is adjustable in height and angle.

SAFARI

3 years
and
above



Inside leg
from 30 to 46 cm

Happy

8 years
and
above




Inside leg
from 50 to 70 cm

Model HAPPY. Cute pink colours in combination with a clear white and a rear tray, ideal to transport toys.

Model SPORTY. Deep dark black topped with fiery red flames and a rear tray, ideal to transport things. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

Sporty

8 years
and
above



Inside leg
from 50 to 70 cm



total
length



total
width



saddle
height



distance
saddle-pedal



step-in
height



height
steering



wheels



total
weight



users
weight

Aqua-Safari	1110	648	500-640	390-520	118	726-838	12" x 2.00	14	40
Happy-Sporty	1294	648	610-760	480-650	133	827-940	16" x 1,75	16,5	60

Sporty - Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

E-bike



17 km/h



36 V / 9 Ah



250 W

Autonomy : depends on the usage,
state of the road, body features
and selected propulsion.

Vintage

VINTAGE - Adults

A classic and stylish retro will get the peoples attention. Brown tyres and a large optional vintage basket completes the picture.

Range of options available at page 128.

Also available in **E-version**, with an integrated electric motor which can be used for propulsion.



Inside leg
from 75 to 108 cm

The steering wheel (different steerings available) is adjustable in height and angle. ■

V brakes at the front with parking lock or an optional disc brake for the rear wheel. ■

Lowered frame for easy mounting/dismounting. ■

Big metal transport tray at the rear. Optional vintage basket. ■

Saddle adjustable in height. ■

LAGOON

LAGOON - Adults

Available in a classic bordeaux colour or a more neutral grey tone. The perfect tricycle for long rides, offering stability and comfort.

Range of options available at page 128.

Also available in **E-version**, with an integrated electric motor which can be used for propulsion.



Inside leg
from 75 to 108 cm



total
length
1808



total
width
768



saddle
height
789-1030



distance
saddle-pedal
639-880



step-in
height
194



height
steering
1121-1309



wheels
24" x 1³/₈



total
weight
27,8



users
weight
120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Tricycles options



01 + B1 - Free pinion



02 - Torpedo brake



03 - 07 - 3 or 7 gears



06 - Fixed pinion



23 - Differential for two rear wheels



B1 - Disc brake rear wheel



BAT 9 - BAT 13 - Battery pack e-Tricycle (9A or 13A)



D1 - Steering lock



E1 - Pushbar (335 mm or 605 mm)



E2 - Front/rear light (batteries included)



E3 - Speedometer (batteries included)



E4 - Vintage basket



E5 - Topcase volume: 30 liters



E6 - Crutchholder



E7 - Can holder (without can)



E8 - Mirror



E9 - Hand straps



E10 - Hand holder with Velcro straps

	width	length
E10 S	60 mm	130 mm
E10 M	90 mm	150 mm
E10 L	120 mm	170 mm




E11 - Closed hand holder



E12 - Comfort saddle PU + plastic



E13-D - Home Trainer



E14 - Front shopping basket



F1 - Pedals with footstrap



F2 - Pedals with feet and toe strap



F3 - Foot support with Velcro



F4 - Foot support + ankle holder



F5 - Adjustable pedal in 2 positions (22 mm)



P1 - Pelvic support



P2 - Back support



P3 - Back and pelvic support (P1+P2)



P4 - Depth adjustable saddle (140 mm)



P5 - Saddle tube with suspension



S1 - Classic steering



S2 - Flat steering



S3 - City steering



S4 - Rectangular steering



S5 - Adjustable steering bracket (-10° à +50°)



S6 - Extended steering

Tricycles options

COMPATIBILITY LIST

		2202	2215	2217	2219	2219 vintage
01 + B1	Free pinion	-	-	•	•	•
02	Torpedo brake	-	•	•	•	•
03	3 gears	-	•	•	•	•
07	7 gears	-	-	•	•	•
06	Fixed pinion	•	•	•	•	•
23	Differential	-	-	-	•	•
B1	Disc brake rear wheel	-	-	•	•	•
BAT 9 - 13	Battery pack e-Tricycle	-	-	-	•	•
D1	Steering lock	•	•	•	•	•
E1	Pushbar	• (S)	• (S)	• (M)	• (M)	• (M)
E2	Front/rear light	-	•	•	•	•
E3	Speedometer	-	•	•	•	•
E4	Vintage basket	-	-	-	-	•
E5	Topcase	-	-	•	•	•
E6	Crutchholder	-	-	•	•	•
E7	Can holder	•	•	•	•	•
E8	Mirror	•	•	•	•	•
E9	Hand straps	•	•	•	•	•
E10	Hand holder with Velcro straps	• (S)	• (S)	• (M)(L)	• (M)(L)	• (M)(L)
E11	Closed hand holder	•	•	•	•	•
E12	Comfort saddle	-	-	•	•	•
E13-D	Home trainer	-	•	•	•	•
E14	Front shopping basket	-	-	•	•	•
E-Tricycle	Motor assistance	-	-	-	•	•
F1	Pedals with footstrap	•	•	•	•	•
F2	Pedals with feet and toe strap	•	•	•	•	•
F3	Foot support with Velcro	• (S)	• (S)	• (M)(L)	• (M)(L)	• (M)(L)
F4	Foot support + ankle holder	•	•	•	•	•
F5	Adjustable pedal in 2 positions	-	•	•	•	•
P1	Pelvic support	• (S)	• (S)	• (M)(L)	• (M)(L)	• (M)(L)
P2	Back support	• (S)	• (S)	• (M)(L)	• (M)(L)	• (M)(L)
P3	Back and pelvic support	• (S)	• (S)	• (M)(L)	• (M)(L)	• (M)(L)
P4	Depth adjustable saddle	• (S)	• (S)	• (M)	• (M)	• (M)
P5	Saddle tube with suspension	-	-	-	-	•
S1	Classic steering	•	•	•	•	-
S2	Flat steering	•	•	•	•	-
S3	City steering	-	-	•	•	standard
S4	Rectangular steering	•	•	•	-	-
S5	Adjustable steering bracket	-	-	•	•	•
S6	Extended steering	•	•	standard	standard	standard

(S) : SMALL - (M) : MEDIUM- (L) : LARGE



VERMEIREN

Opening up new horizons in mobility



easy life

support



Walkers - Rollators	132
Walking sticks	140
Grabbers - Pedalo	142
Commodes	143

Bathing seats	148
Toilet assistance	151
Bars	152

286 E

ROLLATOR

4-wheel foldable rollator. Equipped with a removable basket, plexi table, seat, crutch holder. Parking brakes and height adjustable handles.

Very simple and fast folding system.



780-970 644 605 590 430 11 120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Four Light

ROLLATOR

The four light is a multi-functional 4-wheel rollator with height adjustable handles, a safety belt, parking brakes, a crutch holder and a shopping net. Easy foldable for transport and limited in weight.



870-990 750 550 660 450 8,45 130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Quadri-Light

ROLLATOR

A lightweight, foldable 4-wheel aluminium rollator with parkingbrakes. Ideal for people with (lesser) walking capabilities but who still have the strength and power to use the brakes when they are about to lose their balance.

A smooth plastic seat in combination with a supporting belt when the user wishes to rest.



805-930 700 610 650 290 450 7,25 130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

optional
transportbag
available



276

ROLLATOR



C29
grey

Equipped with a removable basket, a seating cushion, plexiglas tray drawer, anatomic handles, seat and back. Ultra-light foldable model with progressive brakes and parking position.

Provided with extra long pads to rest your elbows.



ECO-light

ROLLATOR



C59
champagne

4-wheel foldable rollator. Equipped with a removable basket, padded seat and crutch holder. A foldable bar provides extra support when seated. Parking brakes and height adjustable handles. Folded and ready for transport in seconds.



1030-1150 750 570 630 440 11 120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



820-950 760 550 600 410 7,6 120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

201 Delta

ROLLATOR

Foldable 3-wheel narrow rollator with a storage bag, removable tablet and shopping basket. Progressive brakes and parking position.



C29
grey



870-990 580 220 630 230 7 130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Quava

ROLLATOR



red

A lightweight, foldable 4-wheel aluminium rollator with parking brakes. Ideal for people with (lesser) walking capabilities but who still have the strength and power to use the brakes when they are about to lose their balance.

A padded seat in combination with a supporting belt when you like to sit down and rest. It's equipped with a crutchholder and the cotton shopping bag can be removed and be carried with you.



805-955 710 540 655 240 450 7,7 136

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Cadeo

ROLLATOR

The Cadeo is a foldable, very light rollator with walk and stop-function. Move forward when walking and push down with your hands to stop. Complete with a removable storage bag and a padded seat with back support.



blue



800 600 540 590 450 5,5 90

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

265 - 266 - 267

WALKING FRAMES

Model 265 demountable

Tubular steel frame painted in blue.
Padded armpit supports which are adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.



C27
blue



Model 266 foldable

Tubular steel frame painted in blue. Padded armpit supports which are adjustable in height. 4 wheels of 75 mm of which 2 with brake.
Foldable for easy transport and storage.



Model 267 demountable

Same features as the mod.265 with the addition of the adjustment of the base height up to 150 mm.



265
266
267

1170-1470
1170-1470
1115-1560

760
760
760

520
520
465-610

580
580
580

min. 300
min. 275
min. 300

13,7
12,5
14

120
120
120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

268 - 269

WALKING FRAMES

Model 268 foldable with table

Tubular steel frame painted in blue.
Padded arm support adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.



C27
blue



Model 269 Foldable with lower arm support

Tubular steel frame painted in blue.
Lower arm support adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.

Options for walking frames

- Table (type 268)
- Lower arm supports (type 269)
- Seat



268
269

830-1130
830-1130

760
760

520
520

580
580

15,5
12,5

120
120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Walkers

250

width: 480 mm
height: 700 - 1000 mm
length: 550 mm
seat height: 620 mm
seat width: 390 mm
weight: 6,5 kg
users weight: 120 kg
colour: champagne



253

width: 560 mm
folded width: 185 mm
height: 795 - 915 mm
length: 670 mm
seat height: 530 mm
seat width: 415 mm
weight: 4,60 kg
removable basket: yes
removable plexi table: yes
users weight: 130 kg
colour: grey



Lyna II

width: 560 mm
length: 470 mm
folded width: 120 mm
height: 780 - 955 mm
weight: 2,65 kg
users weight: 130 kg
colour: inox



Options



Lindsey

width: 590 mm
length: 540 mm
folded width: 120 mm
height: 790 - 970 mm
weight: 2,75 kg
users weight: 130 kg
colour: inox



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.


Streeters

The Vermeiren Streeter was designed to ease the learning and relearning of walking. The locking system on the rear wheels will prevent the Streeter from moving backwards and will encourage the user to move forward.

The system is fully reversible so that your Streeter can be used in both directions. The aluminium frame is fully collapsible which makes it excellent for transport and storage. 4 different colors represent the 4 different sizes.



model MEDIUM




Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.




Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

model EXTRA LARGE



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

model SMALL



model LARGE



options



rubbers



directional wheels

Canes

WALKING STICKS

Foldable models

height: 830 mm - 930 mm
users weight: 100 kg




Beatrix



Bonny



Britney



Zoë tripod

length from 750 mm
to 980 mm
users weight: 130 kg



Mary



Maud



Mona



Monica

Fixed models

height 790 mm - 1020 mm
users weight: 130 kg

Betty



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Canes

WALKING STICKS



Flora 1



Flora 2



Flora 3



Marylin

Two in one. Marylin is a walking stick and an umbrella. Even when it rains, you disassemble the Marylin resulting in an umbrella and a walking stick at the same time.

height from 870 mm to 920 mm
users weight: 130 kg



Electra

Fixed cane with LED lights and loud alert sound
height from 775 mm to 1000 mm
users weight: 130 kg



red



purple



blue



Foldable canes

Flora 1 - Flora 2 - Flora 3

height: 840 mm - 940 mm
users weight: 100 kg

Ashley SMALL & BIG

(adjustable height)
height 1140 mm - 1340 mm (small)
height 1340 mm - 1540 mm (big)
users weight: 120 kg



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Miscellaneous

INDOOR




Nicky

The helping hand Nicky is designed for daily assistance to grab or pick up small objects on the floor.

The front of the helping hand Nicky is equipped with a small magnet. It can help you to pick up small magnetic objects such as keys, coins or cans.

thickness: 12 mm
length: 760 mm
height: 110 mm
weight: 0,15 kg
colors: black, red, green, blue



Twist

Universal rotating cushion, ideal for all types of seats such as your car or your office chair.

diameter: 400 mm
thickness: 50 mm
weight: 1,5 kg
users weight: 120 kg



Kim

Keep moving at home with a hand pedalo. A digital meter keeps tracks of your calories. Fully foldable for easy storage.

width: 400 mm
length: 500 mm
height: 300 mm
disassembled height: 170 mm
weight: 2,3 kg
users weight: 50 kg



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

9060

COMMODOES



C29
grey



550 465 490 470 350 7,2 120



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

9139

COMMODOES



C27
blue



585 450 490 420 450 13,8 100



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



139 E

FEATURES

Easy maintenance

COMMODOES

Commode with wheels and removable arm- and footrests.



C27
blue



560 420 528 420 375 13 120



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

9062 / 9063

COMMODOES



C26
green



600 460 500 470 400 8,7 120
780 (XXL) 600 510 540 400 15 175



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

175 bis

COMMODOES



C26
green



585 470 484 440 360 9,35 95

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Adjustable in height

9063

Height adjustment:
from 410 to 560 mm



9300 series

COMMODOES



C29
grey




C37
dark
blue



Model 9300 big rear wheels

Anatomical shaped back, reclinable and removable middle part of the seat, removable rear wheels 24x1 3/8 PU, directional front wheels 200x50 PU soft, manual brakes on the rear wheels.

Seat width: 430 - 460 - 500 - 530 mm



L06 - removable headrest



L55 - adjustable headrest



B22 - spreader

Accessories for the commode chairs

- Headrest L55
- Headrest L06
- Spreader B22
- Side supports L04
- Wooden table B12
- Plexi table B15

Model 9301 directional rear wheels

Anatomical shaped back, reclinable and removable middle part of the seat, removable front wheels 24x1 3/8 PU, directional rear wheels 200x50 PU soft, manual brakes on the front wheels.

Seat width: 430 - 460 - 500 - 530 mm



9300 series

COMMODOES


Model 9302 wheels of 200x50 mm

Anatomical shaped back, reclinable and removable middle part of the seat, 4 wheels of 200x50 PU with two directionals at the rear, manual brakes on the front wheels.
Seat width: 430 - 460 - 500 - 530 mm



Model 9303 wheels of 100 mm

Anatomical shaped back, reclinable and removable middle part of the seat, 4 wheels of 100 mm with brakes on the rear wheels.
Seat width: 430 - 460 - 500 - 530 mm



9300
9301
9302
9303

660
660
620
560

430-530
430-530
430-530
430-530

260
260
260
260

500
500
500
500

580
580
580
580

450
450
450
450

1070
1070
1070
1070

1100
890
840
850

25,5
26,5
24,5
23,5

100
100
100
100

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Bathing seats

EASY LIFE RANGE



blue

white

All our **Easy life** products are ideal for usage in the shower, bathtub, bathroom or other humid environments.

The special designed rubbers at the end of each leg will create a safe stability for the bench and avoid the sliding of the bench.



Jill



Jinny



Lory

Bathroom chair

Both Lola and Lory are two practical, height adjustable bath stools with rubber feet.

All 4 legs can be adjusted in height by simply pressing a button. For easy storage, all 4 legs can be slid inwards.



Grace

Removable shower seat for the wall.
Easy installation.



Sidney

Bathroom seat with toilet opening

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Bathing seats

EASY LIFE RANGE



blue



white

Lilly

Shower chair with backrest



users
150 KG
weight



users
130 KG
weight

Lola

Bathroom chair with backrest

Both Lola and Lory are two practical, height adjustable bath stools with rubber feet.

All 4 legs can be adjusted in height by simply pressing a button. For easy storage, all 4 legs can be slided inwards.

users
130 KG
weight

Laly

Functional bath- or shower chair. Easy in maintenance and easy to fold after usage.

The legs are adjustable in height depending on your body length.



users
150 KG
weight



blue



white

Melbourne

Bathroom seat with toilet opening



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Bathing seats

EASY LIFE RANGE



white


The transfer benches Kate and Katy facilitate the transfer of a chair to your bath.

Like all products in the range, they are equipped with a non-slip surface. The handle is mounted on the right or left of the bench according to your needs.

Vicky

Transfer board

users
130 KG
weight



Kate

users
150 KG
weight



Vera

Transfer board with 360° rotating disc.

users
130 KG
weight



Magda ALU

The swivel bath tub transfer chair Magda ALU can be adjusted to fit into most bath tubs on the market. Getting in and out of the tub is so much safer with the Magda swivel chair. The seat swivels 90 degrees, so that you can easily transfer in.

users
130 KG
weight



Paris

This high-end bath tub swivel gives u perfect safety and stability. Non-slip surface material and quality tubing for long lasting usage.

users
130 KG
weight



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Toilets

EASY LIFE RANGE



white

Stacy

Commode
Raised toilet seat
Toilet frame
Armrests can be lowered
for easy transfer.

users
130 KG
weight



Tina

Rising toilet seat
Raise the height of the toilet
seat to a comfortable level and
allow the user to easily sit down
and get up from the toilet.

width: 520 mm
depth: 430 mm
height: 380 mm
seat width: 450 mm
adjustable armrests from
5 to 90 mm
height of chair without
armrests: 110 mm
weight: 2,5 kg
users weight: 130 kg



Ina

The bathroom bars are available
in different widths from 300 to 600 mm.
They can be used in almost any
environment.



Liddy

Toilet frame adjustable in height
and ease of use.

width: 560 - 660 mm
depth: 470 mm
height: 630 - 740 mm
weight: 1,45 kg

Installation without drilling.



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

General


EASY LIFE RANGE



white



Izzie L



Izzie


The bathroom bars are available in different lengths from 300 to 800 mm. With exception of the Izzie R and Izzie L, which are bars in L-shape. They are designed as general support elements. They are useable everywhere, including humid environments.

Izzie	30	40	45	60	80 cm
	0,4	0,45	0,50	0,55	0,75 kg

Iris

The foldable toilet bar Iris is designed to help you getting up from the toilet. It can be folded upwards if needed. Additionally, a toiletroll holder is provided.

width: 100 mm
length: 740 mm
height: 230 mm
weight: 1,8 kg
users weight: 130 kg



Izzie 800 mm

Izzie 600 mm

Izzie 450 mm

Izzie 400 mm

Izzie 300 mm

Ivy


Functional and sturdy bathtub rail providing extra support and safety whilst entering or exiting your bathtub. Easy to install and applicable on most types of bathtubs due to its adjustable width (from 60 to 115 mm).




Notes




WHEELCHAIRS




SCOOTERS




RAMPS




CHAIRS



BEDS



PATIENT LIFTS




TRICYCLES



COMMODE CHAIRS - ROLLATORS



GERIATRIC CHAIRS - TABLES




GENERAL SYMBOLS



A form of destructive testing usually performed in order to ensure safe design standards in crashworthiness and crash compatibility for various modes of transportation or related systems and components.

The CE mark, or formerly EC mark, is a mandatory conformity marking for certain products sold within the European Economic Area (EEA) since 1985.

COMFORT AND SHELL CHAIRS



Notes

Handwriting practice lines consisting of 24 horizontal dotted lines.



Notes



Notes

Handwriting practice lines consisting of 24 horizontal dotted lines.



Notes

Handwriting practice lines consisting of 24 horizontal dotted lines.



Notes

Handwriting practice lines consisting of 24 horizontal dotted lines.



Version: 02 2018

All rights reserved, including translation.

No part of this manual may be reproduced in any form what so ever (print, photocopy, microfilm or any other process) without written permission of the publisher, or processed, duplicated or distributed by using electronic systems.

© N.V. Vermeiren N.V. 2018

Vermeiren N.V.
Vermeirenplein 1-15
2920 Kalmthout - Belgium

Tel : +32(0)3 620 20 20
Fax : +32(0)3 666 48 94

info@vermeiren.be
www.vermeiren.com



For more information, contact your local dealer:

www.vermeiren.com
— Quality since 1957 —

Products in this publication are reserved for professional usage.
Subject to technical alterations.