

[illegible]



Version: April 2015

All rights reserved, including translation.

No part of this manual may be reproduced in any form what so ever (print, photocopy, microfilm or any other process) without written permission of the publisher, or processed, duplicated or distributed by using electronic systems.

© N.V. Vermeiren N.V. 2015

Kids

2202 Aqua - Safari

2215 Happy - Sporty

6

6

10

Adolescents

2217 Freedom

14

14

Adults

2219 Lagoon

2219 Vintage

E-version

16

16

18


20

Options

22

How do I choose my tricycle?

The choice of your tricycle is determined by the total length of the inner leg (A).
(total distance, stretched leg, from the folds of the groin to the heel)







Inside leg
from 30 to 46 cm

3 years
and
above






AQUA - unisex model

Decorative figures and bright, shiny colours form the main attraction for young children. The bell is shaped as a dolphin and the rear tray is ideal to transport toys around. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



Range of options available at page 22.

|  |  |  |  |  |  |  |  |  |
|---|---|---|---|---|---|---|--|---|
| total length | total width | saddle height | height saddle-pedal | step-in height | height steering | wheel size | total weight | users weight |
| 1110 | 648 | 498-627 | 390-519 | 118 | 726-838 | 12" x 2.00 | 14 | 40 |

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



SAFARI

2202

Inside leg
from 30 to 46 cm

3 years
and
above










SAFARI - unisex model

Decorative safari figures and bright colours make this tricycle eye-catching to young children. The bell is giraffe shaped and a rear tray is ideal to transport toys. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



Range of options available at page 22.

|  |  |  |  |  |  |  |  |  |
|---|---|---|---|---|---|---|--|---|
| total length | total width | saddle height | height saddle-pedal | step-in height | height steering | wheel size | total weight | users weight |
| 1110 | 648 | 498-627 | 390-519 | 118 | 726-838 | 12" x 2.00 | 14 | 40 |

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Inside leg
from 50 to 70 cm

8 years
and
above

HAPPY - 100% girls










Cute pink colours in combination with a clear white and a rear tray, ideal to transport toys. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



V-brakes

Range of options available at page 22.

|  |  |  |  |  |  |  |  |  |
|---|---|---|---|---|---|---|--|---|
| total length | total width | saddle height | height saddle-pedal | step-in height | height steering | wheel size | total weight | users weight |
| 1294 | 648 | 594-762 | 482-650 | 133 | 827-940 | 16" x 1,75 | 16,5 | 60 |

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.





2215

Inside leg
from 50 to 70 cm

8 years
and
above










SPORTY - Boys will be boys

Deep dark black topped with fiery red flames and a rear tray, ideal to transport things. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



Range of options available at page 22.

|  |  |  |  |  |  |  |  |  |
|---|---|---|---|---|---|---|--|---|
| total length | total width | saddle height | height saddle-pedal | step-in height | height steering | wheel size | total weight | users weight |
| 1294 | 648 | 594-762 | 482-650 | 133 | 827-940 | 16" x 1,75 | 16,5 | 60 |

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.

FREEDOM

WHITE

The Freedom is a sturdy tricycle ideal for demanding adolescents available in black and white. A perfect tricycle for strolling through countryside.

Range of options available at page 22.



total
length

1599



total
width

670



saddle
height

709-950



height
saddle-pedal

586-827



step-in
height

174



height
steering

986-1174



wheel
size

20" x 1,75



total
weight

25,5



users
weight

100

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

FREEDOM

BLACK

- Different options available.
- The steering wheel (different steerings available) is adjustable in height and angle.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- V brakes at the front with parking lock. Optional disc brake for the rear wheels available.
- Big metal transport tray at the rear.



2217

Inside leg
from 68 to 88 cm

14 years
and
above



V-brakes

LAGOON - Adult tricycle

Available in a classic bordeaux colour or a more neutral grey tone. The perfect tricycle for long rides, offering stability and comfort.

Range of options available at page 22.

Also available in **E-version**, with an integrated electric motor which can be used for propulsion.



LAGOON

2219

Inside leg
from 75 to 108 cm

Adults

- The steering wheel (different steerings available) is adjustable in height and angle.
- V brakes at the front or an optional disc brake for the rear wheels.
- Lowered frame for easy mounting/dismounting.
- Saddle adjustable in height.
- Big metal transport tray at the rear.



V brakes



total
length

1808



total
width

768



saddle
height

789-1030



height
saddle-pedal

639-880



step-in
height

194



height
steering

1121-1309



wheel
size

24" x 1.3/8



total
weight

27,8



users
weight

120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.

VINTAGE - Adult tricycle

A classic and stylish retro will get the peoples attention. Brown tyres and a large optional vintage basket completes the picture.

Range of options available at page 22.

Also available in **E-version**, with an integrated electric motor which can be used for propulsion.



Vintage

2219

Inside leg
from 75 to 108 cm

Adults

- The steering wheel (different steerings available) is adjustable in height and angle.
- V brakes at the front with parking lock or an optional disc brake for the rear wheel.
- Lowered frame for easy mounting/dismounting.
- Big metal transport tray at the rear. Optional vintage basket.
- Saddle adjustable in height.



total
length

1808



total
width

768



saddle
height

789-1030



height
saddle-pedal

639-880



step-in
height

194



height
steering

1121-1309



wheel
size

24" x 1^{3/8}



total
weight

27,8



users
weight

120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.



E-version

Electrical model with propulsion support.

The E-version for adults is available for the **LAGOON** and the **VINTAGE** version.
The electric tricycle will allow you to go further, easier and faster.

The operating system is very easy to use. A simple **ON/OFF** switch at the battery and a controller at the steering wheel to get started. It is also equipped with a **BOOSTER** button to give a quick acceleration to the tricycle (for example when you need to start uphill).

The battery pack, located at the back of the tricycle is light and removable, allowing you to easily recharge it at home.

The electrical assistance will power the front wheel of the tricycle.



E-Version

Inside leg
from 75 to 108 cm

Adults



17 km/h



36 V / 9 Ah



250 W

Autonomy : depends on the usage, state of the road,
body features and selected propulsion.

CONFIGURATION

Steering



S1 - Classic steering



S2 - Flat steering



S3 - City steering



S4 - Rectangular steering



S5 - Adjustable steering bracket
(-10° to +50°)



D1 - Steering limiter

Moving



O6 - Fixed wheel



O1 + B1 - Free wheel



B1 - Disc brake for the rear wheel



O2 - Torpedo brake



O3 - 3 gears



E13 - Hometrainer

O7 - 7 gears

Support



P1 - Pelvic support



P2 - Back support



P3 - Pelvic and back support



E9 - Hand straps



E10 - Hands holder with straps



E11 - Closed hands holder

A N D O P T I O N S

Miscellaneous



BAT 9 - BAT 13 - Extra battery E-bike (9A or 13A)



P4 - Adjustable saddle (140 mm)



P5 - Saddle tube with shock absorber



E12 - Comfort saddle PU + plastic



E1 - Pushing handle (335 mm or 605 mm)



E14 - Front basket

Pedals



F1 - Pedals with strap



F2 - Pedals with toes clip & strap



F3 - Foot support with Velcro straps



F4 - Foot support with Velcro straps & ankle support



F5 - Pedals with adjustable length (22 mm)



E2 - Lighting system (batteries included)

Accessories



E3 - Speedometer (batteries included)



E4 - Vintage Basket




E5 - Lockable topcase volume: 30 liter



E6 - Crutch holder



E7 - Drink holder (cannister not included)



E8 - Rear mirror

VERMEIREN N.V.

Vermeirenplein 1-15
B-2920 Kalmthout

Tel: +32 (0)3 620 20 20

Fax: +32 (0)3 666 48 94

website: www.vermeiren.com

e-mail: info@vermeiren.com



Distributed by:

VERMEIREN
Opening up new horizons in mobility