

VERMEIREN



VERMEIREN

**N.V. Vermeiren N.V.**

Vermeirenplein 1 / 15
B-2920 Kalmthout
Tel: +32(0)3 620 20 20
Fax: +32(0)3 666 48 94
website: www.vermeiren.be
e-mail: info@vermeiren.com

L. Vermeiren Ges. mbH

Winetzhammerstraße 10
A-4030 Linz
Tel: +43(0)732 37 13 66
Fax: +43(0)732 37 13 69
website: www.vermeiren.at
e-mail: info@vermeiren.at

Vermeiren Polska Sp. z o.o

ul. Łączna 1
PL-55-100 Trzebnica
Tel: +48(0)71 387 42 00
Fax: +48(0)71 387 05 74
website: www.vermeiren.pl
e-mail: info@vermeiren.pl

Vermeiren Medical Equipm. Co. Ltd.

No. 80, Xingpu Road,
Shengpu Town, S.I.P.
Suzhou - Jiangsu Province
website: www.vermeiren.com.cn
e-mail: info@vermeiren.cn

Vermeiren Deutschland GmbH

Wahlerstraße 12 a
D-40472 Düsseldorf
Tel: +49(0)211 94 27 90
Fax: +49(0)211 65 36 00
website: www.vermeiren.de
e-mail: info@vermeiren.de

Reatime S.R.L.

Viale delle Industrie 5
I-20020 Arese MI
Tel: +39 02 99 77 07
Fax: +39 02 93 58 56 17
website: www.reatime.it
e-mail: info@reatime.it

Vermeiren Nederland B.V.

Domstraat 50
NL-3864 PR Nijkerkerveen
Tel: +31(0)33 2536424
Fax: +31(0)33 2536517
website: www.vermeiren.com
e-mail: info@vermeiren.com

Vermeiren ČR S.R.O.

Sezemická 2757/2 - VGP Park
193 00 Praha 9 - Horní Počernice
Tel: +420 731 653 639
Fax: +420 596 121 976
website: www.vermeiren.cz
e-mail: info@vermeiren.cz

Vermeiren France S.A.

Z. I., 5, Rue d'Ennevelin
F-59710 Avelin
Tel: +33(0)3 28 55 07 98
Fax: +33(0)3 20 90 28 89
website: www.vermeiren.fr
e-mail: info@vermeiren.fr

Vermeiren Suisse S.A.

Hühnerhubelstraße 59
CH-3123 Belp
Tel: +41(0)31 818 40 95
Fax: +41(0)31 818 40 98
website: www.vermeiren.ch
e-mail: info@vermeiren.ch

Vermeiren Iberica, S.L.

Trens Petits, 6. - Pol. Ind. Mas Xirgu.
17005 Girona
Tel: +34 902 48 72 72
Fax: +34 972 40 50 54
website: www.vermeiren.es
e-mail: info@vermeiren.es

Vermeiren India Rehab Pvt. Ltd.

#1489, 29th Main Road, 5th 'c' Cross
BTM 2nd stage, Bangalore-76
Tel: +91 80 41518899
Fax: +91 80 41528899
website: www.vermeiren.in
e-mail: info@vermeiren.in

Index

Manual wheelchairs

Sagitta	6
Jazz	8
Jazz 30°	9
Eclips+	10
Eclips+ 30°	11
Eclips X4	12
750	13
708D	14
V100	15
V200	16
V300	17
28	18
Eclips XXL	19
V300 XL	20
Gemini	21
925	22
708 kids	23
Eclips X4 kids	24
Sagitta kids	25

Comfort

Serenys	28
Inovys II	30
Manual wheelchairs options & gel cushions	32 - 34
Altitude	36
Coraille	37
Version XXL & 917	38
Charme	39

Electronic wheelchairs

Navix	42
Tracer	44
Express	45
Rapido	46
Squod & Squod SU	47
Forest	48
Forest + lift	50
Tracer kids	52
Sunny	53
Forest Narrow kids	54
Springer	55
V-Drive	56
Electronic wheelchair options	57
Ramps	58 - 59

Scooters

Scooters	62 - 72
Scooters options	75

General indoor products

Beds	78 - 85
Patientlifts	86 - 87
Geriatric chairs	88 - 91
Liftchairs	92 - 93
Movement	94 - 99
Toilechairs	100 - 103

Tricycles

2102 - 2103 - 2115 - 2117	106
2119 - options	107

Easy life

Bathing seats & transfers	110 - 111
Toilet accessories	112
Bathroom	113
General (sticks, ...)	114
Walking frames	115

Options

Compatibility tables	116 - 117
----------------------------	-----------





manuals

Manual wheelchairs

Sagitta

Colour



Limited size for transport

Seat and back made from Cordura nylon

Tension-adjustable back and seat

Continuous adjustment of the backrest without extra tools (from -7° to +7°)

Adjustable wheel angles:
0°, 2°, 4°, 6°, 8°
(order only)

V-shaped frame

Pair of anti-tipping wheels - B78



Push handles adjustable in height



Extra cushion



Clothing protectors




Cordura nylon



550-680 570-690 590-710 610-730 630-750 650-770	360 380 400 420 440 460	-7°/+7°	480 (with cushion) 520 (without cushion)	3" 4" 5"	230-280 280-330 330-380 380-430	470-480	320 340 370 390 420 440 460	550 600 650	0° 2° 4° 6° 8°	>9,5	100

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



Clothing protectors
B82S



Clothing protectors
B82



Clothing protectors
B82S + tube



Sport brakes
(scissor brakes)



Standard brakes



Footrest with
aluminium plate (kids)



Footrest with
separated footplates



Monoblock footrest
standard tube




Monoblock footrest
with aluminium plate



Push handles
adjustable in height



Standard
push handles



Wheels 3" translucent



Wheels 4" translucent



Wheels 5" translucent



Wheels 3" black



Wheels 4" black



Wheels 5" black



Spoke protector B85
M3 "Extreme"



Spoke protector B85
M1 "Flowers"



Spoke protector B85
M2 "Flowers"



Spoke protector B85
M3 "Pirate"



Spoke protectors B85
logo Sagitta



Schwalbe Marathon
Plus Evolution



Landcruiser



Schwalbe Right Run



Wheels 24" Off Road



Wheels 24" Spider



Wheels 24" Spinergy



Wheels 24" Viking




Wheels 24" Exclusive



Spoke protectors B85
(transparent)

Jazz

Colour



Light wheelchair with a fixed back

Equipped with Quickrelease rear wheels for easy transport and removable footrests. Completely foldable chair.



Hemiplegic kit



Push version



Lowered version



Amputated version

Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.

590	310	390	210	500	430	430	1070	14		
620	310	420	210	500	430	430	1070	14		
640	310	440	210	500	430	430	1070	14		
660	310	460	210	500	430	430	1070	15		
680	310	480	210	500	430	430	1070	15		
700	310	500	210	500	430	430	1070	16		
						400-510				
										120

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Jazz 30°

Colour



Lightweight wheelchair with a reclining back up to 30° (5 positions)

A perfect straight backrest due to the (removable) push bar. The footrests are removable and foldable. The wheelchair is completely foldable and can easily be transported in a car or other vehicle thanks to the Quickrelease system for the rear wheels.



Push version



Lowered version



Double handrail system



Amputated version

Wide range of accessories and options.


Refer to the accessories page (page 32) and check the compatibility table on page 116.

590	310	390	220	500	480	400-510	430	1070		
620	310	420	220	500	480	400-510	430	1070		
640	310	440	220	500	480	400-510	430	1070		
660	310	460	220	500	480	400-510	430	1070		
680	310	480	220	500	480	400-510	430	1070	21	100

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips+

Colour



Lightweight aluminum wheelchair with a fixed back and reinforced cross

Equipped with Quickrelease rear wheels, foldable and removable armrests.



Amputated version



Hemiplegic kit



Push version

Wide range of accessories and options.


Refer to the accessories page (page 32) and check the compatibility table on page 116.

Front wheel	Front wheel	Front wheel	Front wheel	Front wheel	Front wheel	Front wheel	Front wheel	Front wheel	Front wheel	Front wheel
590	300	390	220	500	420-460	420-440-460	1060	14,5		
620	300	420	220	500	420-460	420-440-460	1060	14,5		
640	300	440	220	500	420-460	420-440-460	1060	15		
660	300	460	220	500	420-460	420-440-460	1060	15		
680	300	480	220	500	420-460	420-440-460	1060	15,4		
700	300	500	220	500	420-460	420-440-460	1060	15,6		
740	300	550	220	500	420-460	420-440-460	1060	15,6		

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips⁺30°

Colour



Lightweight aluminum wheelchair with a reclining backrest of 30°, reinforced cross and Quickrelease wheels

Standard equipped with short armrests for easy table access.



Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.


590	300	390	220	450-520	420-460	900-1040	420-440-460	1000-1170	19	130
620	300	420	220	450-520	420-460	900-1040	420-440-460	1000-1170		
640	300	440	220	450-520	420-460	900-1040	420-440-460	1000-1170		
660	300	460	220	450-520	420-460	900-1040	420-440-460	1000-1170		
680	300	480	220	450-520	420-460	900-1040	420-440-460	1000-1170		
700	300	500	220	450-520	420-460	900-1040	420-440-460	1000-1170		
740	300	550	220	450-520	420-460	900-1040	420-440-460	1000-1170		

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips X4



Colour



Light aluminum wheelchair with adjustable frame

Adjustable in depth, reinforced cross and Quickrelease rear wheels. The footrests can be folded both inward and outward. For optimum comfort, the footrests are adjustable in depth.



Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.

590	300	390	220	390-520	420-460	420-500	1000-1210	14,5		
620	300	420	220	390-520	420-460	420-500	1000-1210	14,5		
640	300	440	220	390-520	420-460	420-500	1000-1210	14,5		
660	300	460	220	390-520	420-460	420-500	1000-1210	14,5		
680	300	480	220	390-520	420-460	420-500	1000-1210	15		
700	300	500	220	390-520	420-460	420-500	1000-1210	15		
740	300	550	220	390-520	420-460	420-500	1000-1210	15		
					400-510					150

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Colour





The 750 provides the basis for a perfect patient positioning. Standard equipped with a backrest adjustment up to 30 degrees and legrests which can be adjusted very smoothly. Transport in a car? Easy, since the 750 is foldable and the arm- and legrests are removable.




Padded flexible back- and seat upholstery

Back adjustable
up to 30°

Lightweight steel
wheelforks



Refer to the accessories page (page 32) and check the compatibility table on page 116.




600	310	390	220	510	480	890	430	1000-1210	21	120
630	310	420	220	510	480	890	430	1000-1210	21	
650	310	440	220	510	480	890	430	1000-1210	21	
670	310	460	220	510	480	890	430	1000-1210	22	
690	310	480	220	510	480	890	430	1000-1210	22	
710	310	500	220	510	480	890	430	1000-1210	22	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

708 D



Colour



Standard wheelchair with a fixed back

The armrests are foldable and removable and equipped with soft, extra long and comfortable padding. The fabric of the seat and back, together with aluminum legrests, will provide optimal comfort while using this chair.



Rear wheels with a Quickrelease system

Sturdy armrests with an easy locking system

Small side wheels on the footplates will protect your walls and furniture while moving around



AMP kit
(amputated version)



708 Delight S
lowered version



T30 kit - push version



Double handrail

Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.

550	280	390	210	510	430	430	430	1100	18	
580	280	420	210	510	430	430	430	1100	18	
600	280	440	210	510	430	430	430	1100	19	
620	280	460	210	510	430	430	430	1100	19	
640	280	480	210	510	430	430	430	1100	19	
660	280	500	210	510	430	430	430	1100	19	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V100



Colour



C58

Metallic blue



C59

Champagne

Model V100

The V100 is a light weight manually propelled wheelchair for temporary or permanent use with a range of settings.

- 4 seat heights
- 2 seat depths
- 9 seat widths
- 3 armrest heights
- 3 positions armrests
- Comfortable seat cushion
- Height adjustable push handles
- Wide variety of options

Also available: version V100D with shopping bag and height adjustable push handles.



Anti-tipping



Optional: plexi table






Brakes



Height adjustable push handles (V100D)

Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.



570

390

220-230-240

470-500-530

420

880-910-940

430-450

1010

starts

130

600

420

220-230-240

470-500-530

420

880-910-940

430-450

1010

from

130

620

440

220-230-240

470-500-530

420

880-910-940

430-450

1010

16,7

130

640

460

220-230-240

470-500-530

420

880-910-940

430-450

1010

1010

130

660

480

220-230-240

470-500-530

420

880-910-940

430-450

1010

1010

130

680

500

220-230-240

470-500-530

420

880-910-940

430-450

1010

1010

130

710

530

220-230-240

470-500-530

420

880-910-940

430-450

1010

1010

175

740

560

220-230-240

470-500-530

420

880-910-940

430-450

1010

1010

175

780

600

220-230-240

470-500-530

420

880-910-940

430-450

1010

1010

175

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V2000



Colour



Model V200

The V200 is a light weight manually propelled wheelchair for temporary or permanent use with optimised settings.

- 5 seat heights
- 2 seat depths
- 9 seat widths
- 3 armrest heights
- 3 positions armpads
- Height adjustable push handles
- Mobile phone pocket under the seat
- Wide variety of options

Also available: version V200D with shopping bag and height adjustable push handles.



Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.

570	390	220-230-240	470-500-530	400	850-880-910-940	440-460	1010	starts	130
600	420	220-230-240	470-500-530	400	850-880-910-940	440-460	1010	from	130
620	440	220-230-240	470-500-530	400	850-880-910-940	440-460	1010	15,2	130
640	460	220-230-240	470-500-530	400	850-880-910-940	440-460	1010		130
660	480	220-230-240	470-500-530	400	850-880-910-940	440-460	1010		130
680	500	220-230-240	470-500-530	400	850-880-910-940	440-460	1010		130
710	530	220-230-240	470-500-530	400	850-880-910-940	440-460	1010		175
740	560	220-230-240	470-500-530	400	850-880-910-940	440-460	1010		175
780	600	220-230-240	470-500-530	400	850-880-910-940	440-460	1010		175

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V30000



Colour



C58

Metallic blue



C59

Champagne

Model V300

The V300 is a light weight manual wheelchair for temporary or permanent use with a maximum on adjustabilities.

- 5 seat heights
- 5 seat depths
- 9 seat widths
- 3 armrest heights
- 3 positions armpads
- 5 adjustable seat angles 0° - 10°
- Tension adjustable back and seat
- Height adjustable push handles
- Mobile phone pocket under the seat
- Wide variety of options

Also available version V300D with shopping bag and height adjustable push handles.



Armrest adjustable in height





Push version



Articulated length auto-adjusting footrest

Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.



570

390

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

starts

130

600

420

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

from

130

620

440

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

15,9

130

640

460

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

1040

130

660

480

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

1040

130

680

500

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

1040

130

710

530

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

1040

175

740

560

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

1040

175

780

600

220-230-240

440-470-500-530

400

850-880-910-940

440-460-480-500-520

1040

1040

175

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

28 Reinforced

Colour



C20
Blue

Model 28

This chair with double cross and fixed back can support a weight up to 150 kg when fitted with an anatomic seat and backrest (110 kg for the standard version).

Equipped with retractable and foldable footrests, desktype armrests that allow easy access to a table and comfortable long armrest padding.

150 KG



Anatomic version (optional)
Max. users weight: 150 kg


Seat and backrest in black nylon

Double cross for additional strength

Small side wheels on the footplates will protect your walls and furniture while moving around

Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.



560	290	390	230	500	420	920	430	1000-1210	18	150
610	290	440	230	500	420	920	430	1000-1210	18	
720	310	500	230	500	420	920	430	1000-1210	19	
790	350	550	230	500	420	920	430	1000-1210	20	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips **XXL**

Colour


C27
Blue

Chairs for wider seatwidths

Aluminum foldable wheelchair with double cross. The Eclips XXL provides a natural comfort with its anatomical seat and backrest. This wheelchair, equipped with reinforced rear wheels, can support a maximum weight of 200 kg.

Easy transfers thanks to the foldable aluminum footplate in one piece.

200 KG

Anatomic seat and backrest



B78
Anti-tipping



B74
Brakes














Reinforced clothing protectors

Foldable, one piece, footplate

Adjustable heel strap

Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.

										
850	375	600	200	530	600	1070	540/590	1140	38,5	200
900	375	650	200	530	600	1070	540/590	1140	39	
950	375	700	200	530	600	1070	540/590	1140	40	
1000	375	750	200	530	600	1070	540/590	1100	40,5	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

V300 XL

Colour



C58

Metallic blue



C59

Champagne blue

Model V300 XL

This chair with double cross can support a maximum weight of 170 kg.
The V300 XL provides a natural comfort with its anatomical seat and backrest.




Also available with double cross:

Model V100 XL and **Model V200 XL**



Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.



	740	530	220-240	440-530	400	840-930	465-505	1085	21	170
V300XXL	770	560	220-240	440-530	400	840-930	465-505	1085	21	170
	810	600	220-240	440-530	400	840-930	465-505	1085	21	170
V200XXL	740	530	220-240	440-530	400	820-1000	490	1050	20,15	170
	770	560	220-240	440-530	400	820-1000	490	1050	20,15	170
	810	600	220-240	440-530	400	820-1000	490	1050	20,15	170
V100XXL	740	530	220-240	470-530	380	890-950	490	1060	21,40	170
	770	560	220-240	470-530	380	890-950	490	1060	21,40	170
	810	600	220-240	470-530	380	890-950	490	1060	21,40	170

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Gemini

Colour



C29
Grey

Model Gemini 1

The Gemini 1 is an aluminum Reha Buggy. The design permits various possibilities of adaptation.

In this standard version we offer a:

- 5-points-security belt
- Spreader
- Headsupport


A buggy which parents adore not only for its excellent functionality, but also for the quality design and above all - the total weight is no higher than 13,5 kg.



Removable rear wheels

Model Gemini 2

Even for larger children we offer a model which is particularly designed towards stability. Compared to the Gemini 1, this model has a steel frame construction without the option of the different seat positions.



Gemini I	630	320	200	500	510-640	250-280-310-340	280-340	770	13,5	35
Gemini II	660	400	200	500	630-750	340-370-400-430	360	980	21	50
Gemini III	740	450	200	500	700-780	390-450	420	980	25	60

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

925

Colour



C8
Yellow



C9
Turquoise
blue

Model 925

The 925 is the wheelchair children love. It is the companion for your child and provides you with the assurance of quality. For example, the robust metal tubing can resist small knocks and special clothing protectors ensure that clothing does not become entangled in the wheels of the wheelchair.

Aluminum
clothing protectors


One-piece footplate

Standard delivered with one-piece foldable footplate adjustable in height and depth.



Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.



480



320




290



120-155



480



360



100-300



330



790



13,4



60

540

320

350

120-155

480

360

100-300

330

790

13,4

60

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

708 kids kids

Colour



C16
Mottled
blue



C29
Grey

Model 708 kids

Compact sturdy wheelchair in steel. Offset of the rear wheels for improved stability. Equipped with the Quickrelease system for fast and easy removal of the rear wheels.



Efficient
brake system

Small wheels on the footplates
protect your walls and furniture
while moving around.



Spoke protectors



Rear wheels with
Quickrelease system



Wide range of accessories and options.

Refer to the accessories page (page 32) and check the compatibility table on page 116.



500
550



310
310




320
370



170
170



510
510



420
420



210-310



350
350



980
980



18
18



70

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Eclips X4 kids

Colour



C34
Yellow



C36
Lilac



C37
Blue

Ultralight kids wheelchair

Constructed in aluminum, the Eclips X4 kids offers you an ultra light wheelchair combined with sturdiness and ease of use. Equipped with clothing-protectors results in a very modern look. Without the wheels (ease to remove) and the footrests, the Eclips X4 kids weighs no more than 9 kilograms making it very easy to transport.

Model T30

T30 Kit - push version



Optional :
Spoke protectors

Padded seat
and back

Quick-release
rear wheels

Aluminum front forks
adjustable in height

Internal and external
collapsible footplates

Armrests that protect
your clothes

Height adjustable
handles (optional)

Wide range of accessories and options.


Refer to the accessories page (page 32) and check the compatibility table on page 116.



500
550



310
310




320
370



140
180-250



390-470
390-470



350-400
350-400



310-330-350-370
330-350-370-390



830-1000
830-1000



17



60

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Sagitta kids

Colour



C53
Orange

C64
White

Active and ultralight kids wheelchairs

Active wheelchair for kids. The Sagitta kids is a mirror image of the Sagitta for adults. And even more important, this wheelchair grows with your child due to the possible adjustments. The frame is constructed in such a way that you can change the seat width.



Sagitta with push bar



Clothing protectors





Footplate adjustable in height



Adjustment of the push bar



Sport brakes (scissor brakes)



460-560
492-592
523-623
555-655

240-340

-7/+7°

370-428
(with cushion)
410-468
(without cushion)

3"
4"

230-430

220-270

280-320

700-750

4°
6°
8°
10°

9

60





comfort

Comfort wheelchairs and
Shell chairs

Serenys

Colour



C20
Blue



C29
Grey



C41
Green





Headrest adjustable in angle,
height and depth

Long backrest: 600 mm (order only)
Short backrest: 500 mm (standard version)

Height ground / seat
adjustable in 3 positions

Easy adjustment
by feet



630
690
750

390
450
510

220-290

420-520
(front)
380-530
(rear)

1020

500/600

19°/113°

-5°/50°

-2°/15°

450-500

1170

120

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

Model Serenys

Very comfortable chair for prolonged seating. Equipped with memory foam cushions, wraparound armrests, a deformable headrest and legrests with retractable cushions facilitating transfers. Smooth tilting back and forth due to two springs. Possibility to lower the seating to 450 mm (with standard cushion) without changing the wheel size.



The legrests are fully adjustable and the footrests with inclination will give an ideal angle for the feet at any given time.



Comfortable, height adjustable armrests. Retractable for easy transfers.



Double gas cylinders providing assistance when the patient wants to get up.



This chair is the most comfortable one on the market for prolonged seating due to its memory-foam seat cushion. Foam density: 41 kg/m³.

Seat angle and seat height are adjustable



The legrests are adjustable in height.



The calf supports are adjustable in height, width and depth.



The calf supports are collapsible to the back to ease transfers.

Inovys II

Colour



C55
Inox

Model Inovys II

Very comfortable wheelchair with multiple settings.

- Adjustment of the armrests in length and height
- Adjustment of the seating depth
- Adjustment of the footrests in width
- Possible to lower the chair without changing its original wheels



Easy to disassemble




Two foot pedals



Side support



Table



630
680
730

400
450
500

220-320

460-530

1020-1200

550-585

10°/120°

0°/54°

0°/21°

460-520

1010-1050

135

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible



Crash Test Approved

Model **Inovys II Maxi comfort**

Standard equipped with:

- Wide and gutted armrests
- Comfortable headrest L58C
- Comfortable legrests BZ8C



Optional comfort headrest



Transit kit with 16" wheels



Removable armrests.
Adjustable in height and width.
Seat adjustable in depth.



Adjustable hemiplegic
armrest



Headrest
adjustable in 3D



Options



B02
long armrest



B03
short armrest



B05
armrest adjustable
in height



B06
standard footrest



BA6
shortened footrest



B08/07
adjustable footrest
08 : short - 07 : long



BK6
footrest in one
piece



BK7/BL7
AMP legrest



BB6
adjustable short
footrest



BZ10
legrest with
side corners



B12
wooden table



B14
backrest extension



B17
toilet seat




B18
heel support



B19
foot straps



B20
safety belt




B21
wide calf support XL



B22
spreader



B23
cushion: 40 mm



B26
cushion: 20 to 60 mm



BZ7/BZ8
articulated length
auto-adjusting footrest



B30
handrail for
arthritis



B31
crutch holder



B40
brake extenders



B42
shopping net



B80
Quickrelease
system



T30
pushing
wheelchair kit



B52
serum holder



B85
spoke protectors



B66/o
gutted rotating
armrest



B73
arthrodesis seat



B74
brakes



B78
anti-tipping



B79
document holder



B90
push bar



B91
anti-theft slat



L44
tension adjustable
backrest




L55
shapeable headrest



L04
side support



B86
travel wheels




B115
plexi table



Options



L14
anatomic back
(soft)



L15
anatomic seat
(soft)



L14
anatomic back
(skai)



L15
anatomic seat
(skai)




L24
preshaped back
(skai)



L25
preshaped seat
(skai)



L34
preshaped back
(soft)



L35
preshaped seat
(soft)



BL6
one-piece footrest
shortened



BN6
shortened footrest




BZ9
one-piece adjustable
legrest



BZ10
BZ9 legrest with side
corners of 70 mm



B41
pump



B37
60 mm cushion with
oval opening



B58
safety belt



B84
plastic spoke
wheels



B101
anti-theft system
with money-coin



B112
hemiplegic armrest



AMP
AMP kit: offset of
the rear wheels



HEM2
hemiplegic kit:
double rims



S1
lowered
seating kit



L02
anatomic seat with
elevating system in front



L13
backrest hooks



B83
cannister holder

Essentials



V10

Gel protective cover
60 x 240 mm (pair)

Essential gel cushions

There are many cases where there is a need for additional protection of certain body parts while using a wheelchair. The Vermeiren company is proud to offer a wide range of different gel comfort pads that are easy to apply and can be used in almost any situation.



V12 - V14

Protective gel covers for
tubing 6" - 8" (pair)
120 x 160 mm - 120 x 200 mm



V16

Gel heel protectors (pair)
120 x 140 mm



V18

Velcro gel pad for walker or crutch
(pair) 120 x 130 mm



V20

Gel protection for a
safety belt (2 parts)
325 x 80 mm and 230 x 70 mm



V22

Gel protection for gutted
armrest (pair) 110 x 340 mm



V24

Gel leg side protection
(pair) 310 x 90 mm



V26

Gel strap support for legs
2 pieces 400 x 190 mm



V28

Gel plate footrest
(pair) 155 x 210 mm



V30

Adhesive oval pad
(pair) 80 x 100 mm

Essentials

**V32**

Adhesive oval pad
(pair) 100 x 200 mm

**V34**

Adhesive oval pad
(pair) 120 x 300 mm

**V36**

Adhesive ovoid pad
120 x 70 mm

**V38**

Gel armrest pad
(pair) 90 x 300 mm

**V40**

Gel protective knee pad
(pair) 100 x 100 mm



Altitude


Colour






Model Altitude

The Altitude is basically equipped with an 80 mm thick seating cushion (backheight: 770 mm). The height of the back and armrests are depending on the thickness of the seating cushion (60-80-100 mm).

Optional anti-decubitus cushion available.



Adjustable seat height
(from 450 to 560 mm)



670

380
400
420
440
480
550
600

230 (cushion 60mm)
250 (cushion 80mm)
270 (cushion 100mm)

450.470.490
510.520.530
540.550.560

1190
1215
1240

750
770
790

-4/36°

460

1120

36

120

790 (XXL)
790 (XXL)

230

540.550.560

1240

780

500

1120

37

150

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Coraille

Colour



Model Coraille


Shell chair for indoor / outdoor usage with 4 wheels: 2 big rear wheels and 2 front wheels.

Optional anti-decubitus cushion available.



Optional table


Leg rest extension
with automatic
self-elongation



Footrest plate mounted
on a spring



Adjustable seat height
(from 450 to 560 mm)



670	380 400 420 440 480 550 600	230 (cushion 6cm) 250 (cushion 8cm) 270 (cushion 10cm)	450.470.490 510.520.530 540.550.560	1190 1215 1240	750 770 790	-4/36°	460	1120	36	120
790 (XXL) 790 (XXL)	480 550 600	230	540.550.560	1240	780		500	1120	37	150

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


Options for Altitude/Coraille



B12 - Removable and foldable table with telescopic adjustment
(seat width 380 mm to 480 mm)



B52 - Serum holder



Optional separate footrests
(request during order)



Foot brake at the bottom

Also available: Optional anti-decubitus Dartex cushion (80 mm)
(available in green, orange or purple)

XXL Altitude Coraille

Colour



dartex green



dartex orange



dartex prune



skai light purple

Version XXL

Shell chair Altitude and Coraille XXL
Seat width 550 mm or 600 mm for an overall width of 790 mm
Maximum users weight: 150 kg

Only available with separated footrests.

XXL
MAX. USERS WEIGHT : 150 KG



917 Shell chair frame



917 T30

Vermeiren delivers also a shell seat frame which is completely adaptable (model 917). The frame is equipped with a mounting system: the seat can be mounted or removed from the frame without any difficulty.

The armrests are adjustable in height, width and angle. In addition, you can choose between classic legrests or a legrest in one piece. The adjustment of the backrest by jacks is optional.

Charme

Colour



Model Charme

The Charme shell chair is the result of design, comfort and aesthetics for people with a weakened muscle tone. The chair allows the user to easily position himself in a resting position and keep his autonomy. It integrates in the patients home due to its looks and design and adds a therapeutic value. Coordinated tilt of seat and legrest for the perfect positioning without complicated adjustments.



Charme Deluxe (black skai)



Comfortable removable headrest

Ergonomically shaped legrests for comfort (automatic extension)


Safety harness

Seat height adjustable in 2 positions

4 swivel casters for effortless movement



Electrical version



740	380 400 420 440 480	230 - 270	470 - 495	1220 - 1245	750 - 790	-10/30°	480	740	25	120
-----	---------------------------------	-----------	-----------	-------------	-----------	---------	-----	-----	----	-----

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible





electro

Electronic wheelchairs

300 mm lift for a better social integration

The lift allows me to interact with people around me, enhancing my mobility and social integration.

300 mm lift for a better accessibility

- Kitchen cupboards
- Cash dispensers
- Bar
- etc...

Extra option: LIFT



Navix

Colour



C30
Gris
metallic



C68
Clear blue
metallic



C69
Brown
metallic



Better maneuvering with a compact front-wheel drive

The Navix allows for sharper turning indoor and for better overcoming obstacles like thresholds when entering the local grocery store.

The user friendly Navix is adjustable to individual needs and allows for easy transfers to bed, bath, toilet, car etc...

The compact Navix fits in a small car for even more mobility.


- Better overcoming obstacles
- Better taking sharp turns
- Easy to reassemble
- Total width for Navix limited to 570 mm
- Ground clearance of 60 mm
- Adjustable to your needs



Adjustable seat and back



Adjustable armrests in height and width



Total width indication






Better grips on outdoor surfaces with profiled large front wheels



2x12V 38Ah	2x220W	6 km/h	20 km	9°	60 mm	1250 mm
---------------	--------	--------	-------	----	-------	---------

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



570
570
620

400
450
500

260-350
260-350
260-350

400/450/500
400/450/500
400/450/500

380-500
380-500
380-500

580
580
580

332-402

945
945
945

79 kg
79 kg
79 kg

130
130
130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Tracer

Colour



Electronic and made-to-measure

The control unit with microprocessor gives an extra dimension in using an electronic wheelchair: it is easy to control, with a minimum of effort and maximum safety.

The motors in the Tracer are the latest generation of motors. The 200/300 Watts enables you to go up slight inclines without problems. The rechargeable and maintenance-free batteries (2x12V/70 Ah) provide enough power for about 45 km.



Anatomical seat and back



Optional:
Koenig version



DX2 controller with
color LCD display




Optional: seat and back
in Soft version



2x12V-70Ah 2x200W 6 km/h
10 km/h 35 km 10°/17% 70 mm 1400 mm

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



670 390 215-265 550 380 490 360-510 1020 129 kg 125
670 450 570 570 450 530 1220 137 kg 125
750 550 570 450 530 1220 137 kg 150

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Express

Colour



C20

Blue

C29

Blanc

C29

Grey

Solid and effective


Electronic wheelchair with a double cross, equipped with footrests and armrests adjustable in height and width.

Front wheels 200 x 50 and rear wheels 12^{1/2} x 2^{1/4}.
Powerful engine of 200W and 38A/h batteries.



A-Series controller

- optimal performance
- low power consumption
- ergonomic



Very compact chair when folded, which facilitates the transport.



2x12V-38Ah 2x200W 6 km/h 30 km 10°/17% 15 mm 1400 mm

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



560	380	390	min. 200	530	460	450	400-510	1050	62	90
580	380	420	max. 280	530	460	450		1050	62	
600	380	440		530	460	450		1050	62	
620	380	460		530	460	450		1050	62	
640	380	500		530	460	450		1050	62	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Rapido

Colour




Model Rapido

Programmable electronic wheelchair for indoor/outdoor usage.

Tension adjustable back

Colour



FREE Option

Seat and back in Nylon Cordura

model AA2

FREE Option

Lighting

Optional: manual adjustment of the seat angle



Folded for transport



2x12V-38Ah 2x200W 6 km/h 30 km 10°/17% 40 mm 1400 mm

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



590	360	410-450	min. 200	550	440	500	400-510	1050	81	120
640	360	460-500	max. 250	550	440	500		1050	81	120
700	360	520-570		550	440	500		1050	81	120

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Squod S.U.

Colour

Model Squod S.U.

C55
Inox

Electronic stand-up wheelchair



Upholstered seat
and back

Model Squod

Colour

C30 C41
Diamant Green
grey



Shark 2
Controller

Foldable
footrests

In- and outwards
foldable legrests

Comfort

Ideal compromise for indoor / outdoor usage. The standing function will adapt to each situation.

The Squod is equipped with standard lighting.



Adjustable back 30°



Armrests adjustable in
height and width






Anti-tipping wheels



2x12V-70Ah 2x200W 6 km/h 30 km 6°/10,5% 70 mm 1400 mm
2x12V-50Ah width 39

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.

**Squod**590
630
700390-430
440-490
500-550

200-250

530

430

510

1070

111

130

Squod SU

650

450

min. 220
&
max. 310

580

min. 400
&
max. 450

510 / 610

1250

121

130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible



A powerful engine of 2 X 350 W for the 10 km/h version and 2 X 200 W for the 6 km/h version. Combined with maintenance free AGM batteries with a capacity of 70 A/h.

Stability axle



Seat depth adjustable on the frame of the wheelchair.



Forest



Colour



C21 Red C53 Orange C55 Inox RAL9207 White

Ideal for indoor and outdoor usage.

The FOREST is ideal for crossing obstacles in front of you. The 4 wheels of the Forest will **always be in contact with the ground** in spite of uneven terrain up to 120 mm in height differences.

New ventilated coating for the seat and the back avoiding any problem of perspiration.

Model Forest AA2



Adjustment of the seat and back by electric jacks, headrest and BZ7 manual footrests. (BZ7E is the electric version)



Equipped with the Dynamic Shark 2 electronic controller



Programmable DX2 controller available



Shock absorbers and anti-tipping at the back



Optional fenders


Armrests adjustable in height



2x12V-70Ah 2x200W or 2x350W 6 km/h or 10 km/h 35 km 10°/17% 70 mm 1400 mm

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



580 390-440 200-250 550 430-490 550 360-510 1170 127 kg 130
650 450-490 500-600 1200 127 kg 130
650 1200 127 kg 130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

- **Adjustable to your needs**
- **Excellent outdoor and indoor performance**
- **Stability and easy to handle**
- **Quality & reliability**
- **Attractive design**

300 mm lift for a better ACCESSIBILITY

- Kitchen cupboards
- Cash dispensers
- Bar
- etc...



Forest Lift

Colour



300 mm lift for a better social integration

The new Forest Lift allows me to interact with people around me, enhancing my mobility and social integration and is adjustable to your needs. Attractive looks and quality design. Robust quality of parts and frame supports efficient use in demanding outdoor and indoor situations. Excellent indoor and outdoor performance for more independence and mobility.



- Excellent maneuverability
Motor 2 x 200W or 2 x 350W
- Excellent outdoor stability. The stability axle keeps the wheels in touch with the ground on uneven terrains up to 120 mm.
- Driving range: 30 km

Armrests adjustable
in height and width

Reclinable seat
and back

300 mm

Legrests adjustable
in length

Footrests adjustable
in length



Better social integration



Legrests adjustable
in width



Footrests adjustable
in length




Adjustable seat depth



2x12V-70Ah 2x200W or 2x350W 6 km/h or 10 km/h 30 km 10°/17% 70 mm 1400 mm

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



580 390-490 200-250 530-830 420-480 510-610 4°-31° 830 155 kg 130
650 440-540 (without seatcushion) 500-600 (without footrests)

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Tracer kids

Colour



C34
Diamond
yellow

Programmable

Mostly outdoor usage wheelchair with electronic adjustable seat and back, equipped with shock absorbers for maximum comfort. Standard lighting.



Comfort

Shaped seat and back for ideal comfort. Armrests adjustable in height and width.




DX2 controller with a large high definition LCD display and 5 customizable programs



2x12V-70Ah 2x300W 6 km/h 35 km 10°/17% 70 1400

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.




640 320 350 370 170-210 550 350-390 500 270-380 1020 130 kg 125

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Sunny

Colour



C4

Red

C8

Yellow

Electronic and made-to-measure

The sunny is an electronic wheelchair for children which has a microprocessor. It is easy to operate, with a minimum of effort and maximum safety. Due to the 5 programs, the wheelchair adapts itself to the specific needs of the child.



Adjustable position for the controller

Adjustable armrests in height and depth

Footplate adjustable in height



Adjustable seat depth



Armrests adjustable in height



2x12V-38Ah 2x200W 6 km/h 30 km 10°/17% 15 mm 1400 mm

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



490	400	320	225-150	530	300-370	490	100-300	800	75	60
510	400	350	225-150	530	300-370	490	100-300	800	75	60
540	400	370	225-150	530	300-370	490	100-300	800	75	60

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Forest Narrow kids

Colour



C55

Inox

C44

Orange

C57

Red

Excels on uneven surfaces

The Forest Kids is an electronic wheelchair for children, used for indoor and outdoor and ideal in overcoming obstacles. The 4 wheels are always in contact with the ground even on rough surfaces up to 120 mm.



Electronic tilting of seat and back.

Height adjustable
pushing handles

Armrests adjustable
in height

Shark 2 Controller

Legrests adjustable
in length

LED lights in front
and back



Legrests adjustable
in width



Stability axle



Excels on uneven
surfaces



2x12V-70Ah 2x200W

6 km/h

30 km


10°/17%

70

1400

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



580



370




170-230



480-520



340-400



400



830-1070



85




130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Springer

Colour



C67 Orange C66 Blue

- Easy outdoor & indoor maneuvering
- Seat lift for better access & social interaction
- Great looking & highly configurable



Seat LIFT up to 350 mm for

- Better social interaction at similar level with friends & family at home, at school, ...

Better access to places difficult to reach

- Shelves at school, in library, at home, in telecom shop, ...
- Cupboard, drawers, ... in kitchen, living room, bath & bed room, ...

Anti-tipping wheels at the back

LED lights in front and back

4 shock absorbers

Profiled wheels for a better grip on uneven terrains

Ground clearance 80 mm



2x12V-38Ah

2x200W or 2x350W

10 km/h

23 km


10°/17°

70

1400

Wide range of accessories and options.

Refer to the accessories page (page 57) and check the compatibility table on page 116.



600

260-380

160-270

530-810*
*with lift

310-370

370-490

15°

-5/+55°

1010

80
90 with lift

100
50 with lift

Dimensions en mm, poids en kg. Tolérances dimensionnelles et pondérales admises : 15 mm, 1,5 kg et 1,5°

V-drive

Comfort

The VDrive system simplifies your movements with a manual wheelchair. The V-Drive takes over all the work from your pushing assistant. Taking a steep hill or inclinations, bumps in the roads or just a small journey outside... all this is just ideal for the V-Drive.

Standard version
for people up to 135 kg
(weight chair included)

H.D. version
for people up to 200 kg
(weight chair included)



Montage

A very simple installation is available for left or rightside mounting of the V-Drive.

Disassembly of your V-Drive for transport inside your car can be done in a minute.

Disconnect the battery pack and the driving system and you are ready to go.

Delivered in a transport bag



standard
H.D.

390-480
500-600

6 km/h
6 km/h

8°
8°

15 km
15 km

2 x 12 V 12 Ah
2 x 12 V 12 Ah

150 W
200 W

(without batteries)

11,8
8,2

135
200

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Options

Options and controllers

Vermeiren offers a wide choice of options and controllers which are adaptable to your electronic wheelchair. Everything is designed to improve comfort and facilitate the usage of the patient and companion. Vermeiren solutions allow you to customize your wheelchair so that everyone can find the solution that will push the limits of mobility.



Vermeiren may at any time customize your wheelchair with adaptations whose only limit is your imagination. As such, each user is entitled to access new technologies and push the limits of mobility.

Controllers



SE84
Board controller



SE79
Hand controller



SE56
Finger controller



SE64
Remote controller



SE82
Button controllers



Reinforced controller



Mini joystick



SE09
Double controller



SE04
Chin controller



SE04e full
Chin controller (electrically retractable)



SE02
Table controller

Options



B12
Plexi table



BK7/BL7
Legrest for amputated



BZ7/BZ8
Legrests 7: long - 8: short



BZ7/BZ8
Electrical legrests



L04
Side support



L55
Shapeable headrest



SE35
Pavement climber



SE52
Adjustable back 30°



AA2
Electrical adjustable back



SE39
Electrical seat and back



SE63
Softball joystick



SE42
Lift



SE67
Seat and back KOENING version



Joystick adaptations




Ramp Kits systems

Model n°4



Model n°1



Model n°2



Model n°3







Kits n°2, 3 and 4 are delivered with an installation CD.



Determine the height to cross and then select the proper ramp kit.

From 10 to 150 mm, this modular system allows you to fine tune the kit to your desired height and support up to 200 kg.

				
				kg
Kit 1	260	750	10 - 40	1,3
Kit 2	510	750	45 - 75	4,3
Kit 3	760	750	75 - 115	9,2
Kit 4	1010	750	115 - 150	15,8

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


691 Ramps

Small model

Large model





Telescopic aluminum ramps with non-slip surface for maximum safety.



Customized Vermeiren products

Vermeiren offers you to create your own customized ramp system.
Kindly contact our sales department..



691
691

1060 - 2000
1660 - 3200

190
190

50
50

7 x 2
10 x 2

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible





electro Scooters

Antares 3 - Antares 4

Colour



Blue

Model Antares 3 / Antares 4

Scooter for mixed usage. Comfortable and very practical, the Antares can be disassembled into 5 parts in seconds without additional tools.



Anti-tipping wheels

Model Antares 3 – scooter with 3 wheels

Shopping basket



Swivel seat

Battery

Options and accessories page 75.

Model Antares 4



Antares 3
Antares 4

6

6

10 - 12

10 - 12

> 8°


> 8°

54

54

1475

2460



480
480

980
980

965
970

200 x 50 / 3
200 x 50 / 4

500-540
500-540

30
31
(without batteries)

115
115

2 x 12V / 12 Ah
2 x 12V / 12 Ah

270 W
270 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Venus 3-Venus 4

Colour

Model **Venus 4**



Easy disassembly
for transport

Swivel seat

Colour

Model **Venus 3**



Swivel seat

Adjustable steering unit






Shopping
basket

Swivel seat

Colour












Options and accessories page 75.

					
(3)	6	16	6°/10%	50	1330
(4)	6	16	6°/10%	50	2420
(4s)	6	12	6°/10%	55	2470

Model **Venus 4 Sport**

4-wheel scooter for mixed usage. Comfortable and very practical, the Venus can be disassembled into 4 parts in seconds without additional tools.

									
(3)	490	910	975	200 x 50 / 3	440-540	49	115	12V / 12 Ah	336 W
(4)	490	915	1030	200 x 50 / 4	440-540	50	115	2 x 12V / 12 Ah	336 W
(4s)	550	900	1040	220 x 80 PU/4	440-540	52	115	2 x 12V / 14 Ah	350 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Carpo 3 - Carpo 4

Colour



Model Carpo 3 – scooter with 3 wheels

Model Carpo 4



Powerful outdoor comfortable 4-wheel scooter.
Electronic programmable. A reliable partner for your trips.



Options and accessories page 75.

**Carpo
XD model***

* Equipped with AGM batteries (2 x 12V 70Ah)
and a more powerful engine for hilly regions.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Carpo 4 Limited edition

Colour



Brown

Model Carpo 4 Limited edition

Powerful and comfortable 4-wheel scooter for outdoor usage. Electronic programmable. The Carpo 4 Limited edition is standard equipped with a LCD display with odometer. This scooter gives you the ultimate driving pleasure.

Comfortable adjustable seat with headrest and pouch



Rear mirror

Extra strong shopping basket

Front and rear lighting and indicators

Bumper

Anti-slip floor mat


Options and accessories page 75.



17



45



9°/15%




150



2900



660



1200



1450




3.00-8



450-680




130



150



2 x 12V / 70 Ah



750 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible



Carpoe

Champagne

Model Carpo Deluxe

Powerful and comfortable 4-wheel scooter for outdoor usage. Electronic programmable. A total new design giving you that very futuristic look.

Comfortable adjustable seat with headrest



Rear mirror



Modern designed steering column

Front and rear lighting and indicators

Bumper

Anti-slip footbase

Options and accessories page 75.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Ceres 3 - Ceres 4

Colour



Black



Blue

Model Ceres 3 – scooter with 3 wheels

Model Ceres 4

Carpo's little brother, the Ceres 4 is a comfortable scooter for your daily usage.



Options and accessories page 75.

km/h	km			
C3 9	30	8°/14%	100	3020
C4 9	30	9°/15%	100	2750

				kg	kg	kg	kg	kg	kg
C3 610	1160	1270	4.00-5	420-690	92	135	2 x 12V / 50 Ah	470 W	
C4 610	1160	1310	330 x 100	420-690	115	135	2 x 12V / 50 Ah	470 W	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Sportrider

Colour



Black

Model Sportrider

The Sportrider is a powerful three-wheeled scooter with retro design, motorcycle style front fork and wheels with chrome mudguards.

Height adjustable headrest

Rear mirror

Hand brake

XL shopping basket

Large headlight

Fenders

Anti-slip rubber floor mat

Options and accessories page 75.



15



50



9°/15%




100



1573



780



1120



1650




2.75-14



490



188,5



182



2 x 12V / 100 Ah



1300 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Triton

Colour





Blue

Model Triton

Powerful outdoor 4-wheel scooter with an extra wide seat.
Maximum users weight: 200 kg



Options and accessories page 75.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Oberon

Colour



Red

Model Oberon

Double seated 4-wheel scooter

Adjustable headrest

Shopping basket

Double seat




Options and accessories page 75.



6



30



8°/14%




150



3500



920



1090



1370




330x100 air/4



590




125



200



2 x 12V / 35 Ah



2 x 520 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Mercurius 3 - 4

Colour



Champagne

Model Mercurius 3








Adjustable seat










Model Mercurius 4



Shopping basket

Options and accessories page 75.

					
	km/h	km	°/15%	cm	kg
M3	15	45	9°/15%	60	1800
M4	15	45	9°/15%	100	3200

									
	mm	mm	mm	front rear	mm	kg	kg	V Ah	W
M3	700	1300	1310	4.00-5 air/1 - 3.30-8 air/2	450-525	142	150	2 x 12V / 70 Ah	1490 W
M4	680	1300	1450	3.30-8 air/4	450-525	143	150	2 x 12V / 70 Ah	1490 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Cabine

Colour



Red



Model **Oberon**

4-wheel all-weather-scooter.

Jupiter 4 fast






Colour



Blue



Options and accessories page 75.

					
CAB	12	45	9°/15%	100	3200
J4 F	15	30	8°/14%	100	2740

								
CAB 850	1620	1500	300x100 air/4	450-525	173	150	2 x 12V / 70 Ah	1490 W
J4 F 620	1300	1280	2.50-8 air/4	480-570	112	120	2 x 12V / 50 Ah	1490 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Saturnus 3-4

Colour



Blue

Model Saturnus 3



Adjustable headrest

Model Saturnus 4

















Shopping basket

Bumper

Adjustable seat

Options and accessories page 75.


					
	km/h	km			
S3	6	25	6°/10%	40	1910
S4	6	25	8°/14%	40	2900

									
	mm	mm	mm			kg	kg		W
S3	550	1020	1130	front 230x65 air/1 - rear 260x85 air/2	360-530	74	115	2 x 12V / 36 Ah	560 W
S4	580	1150	1350	330x100 air/4	360-530	82	120	2 x 12V / 36 Ah	560 W

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Options

Vermeiren offers a wide selection of accessories for your scooter.



Crutch holder

01



Rollator holder

02



Raincover

03



Rear mirror

04



Protective cover

05



Cannister holder

06



Luggage box

07

Attention : Options 1 are 6 incompatible.

option	Antares	Cabine	Carpo	Ceres	Jupiter	Mercurius	Oberon	Sarturnus	Sportrider	Triton	Venus	Venus 4 sport
01	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
02			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>			<input type="checkbox"/>
03			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
04	<input type="checkbox"/>	S	S	S	<input type="checkbox"/>	S	<input type="checkbox"/>	<input type="checkbox"/>	S	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
05			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		
06	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
07		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>

S standard

☐ compatible option





Beds

Indoor care equipment

Luna

Model Luna

The Luna bed is a stylish, high quality 4 section profiling electric bed with a user handset for convenience and full length wooden side rails for safety. The bed comes complete with a lifting pole for extra user movement and independence. The multi adjustable 4 section mattress platform provides comfort and assists correct posture.



Delivered with a transportkit




Highest position



Anti-trendelenburg



Delivered mounted on a transportkit.
Within minutes you can install your bed wherever you want.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Luna Basic Basic

Model Luna basic

An economy three motor version of the Luna bed with similar functions apart from the knee buckle facility and covered bed ends. Legs can be raised manually or, with the four motor version, electrically. Easy to use in all homecare environments.




Delivered with a transportkit



Anti-trendelenburg



Delivered mounted on a transportkit.
Within minutes you can install your bed wherever you want.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


Luna Deluxe

Model Luna Deluxe

A combination of manufacturing excellence, aesthetics and functionality are offered in the optimal version of the Luna bed. Contemporary bed ends together with the top rated functions of the Luna bed family provide an elegant, yet robust, addition to the range.



Delivered with a transportkit




Anti-trendelenburg



Adjustable self-raising device

Delivered mounted on a transportkit.
Within minutes you can install your bed wherever you want.




Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


Luna X-low X-low

Model Luna X-low


The Luna Ultra Low is a 4 section electric profiling bed that is ideal for those confused or disorientated patients who may be prone to falling out of bed. With the variable height position in its lowest setting and with a soft foam roll-out floor pad in place the dangers of injury from unscheduled bed exit can be minimised.




Delivered with a transportkit




Optional: bedsheet arch



Anti-trendelenburg



Delivered mounted on a transportkit.
Within minutes you can install your bed wherever you want.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Pratic

Model Pratic

The Pratic bed is equipped with a head-rest with electronic controls. A simple touch on a keypad enables you to raise it for reading, watching television or to take your meals. The legrest can also be adjusted to any desired position.

Options

- Centralized brake
- Louis-Philippe pannels
- Independent serum holder
- ...



Optional Louis-Philippe pannels
light oak or dark oak.




Delivered with a
transportkit




Wooden side barriers

Adjustable selfraising device
with serum holder.
Possible to install on the left or
right side of the bed.

Next to the side rails is an
additional rail to keep the
mattress into place.



Delivered mounted on a transportkit.
Within minutes you can install your bed wherever you want.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Arbo

Model Arbo

The Arbo bed is not only extremely functional but fulfils all known technical demands. We use the Linak motor to guarantee reliability and an optimal life expectancy. The whole production process is done in-house thus guaranteeing a great service at a very competitive price.

Options

- Centralized brake
- Independent serum holder
- ...



Delivered with a transportkit




Adjustable backrest with auto regression



Metal siderails version



Delivered mounted on a transportkit.
Within minutes you can install your bed wherever you want.



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Interval


Model Interval

The Interval bed has an electrically adjustable head support. By simply pressing a button, the head support moves into a more vertical position, enabling you to read, eat or watch television. The footrest can be placed in the correct position manually.


This high/low bed has 4 wheels, each with a separate brake, therefore guaranteeing stability. The strong and flexible coated baseboard is easy to maintain. The bed can be assembled and disassembled quickly and without tools.

Options

- Bed with electrical knee buckle
- Table
- Bedside table
- Cupboard - table - chair
- Cover to transport the wooden panels
- Wooden base frame
- Centralized brake



Centralized brake



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Interval XXL

Model Interval XXL

The 3-functions Interval XXL bed has electrical adjustments and is available in 1200 and 1400 mm. frame with a grid mattress base. provided with an electrical

whereas the footpart can be

Delivered in a white
The legpart is
adjustable knee-buckle
adjusted manually.

Max. users weight : 270 kg.



Metal siderails version



1200x2000	1330	2080	370	805	-	30°	70°	30°	240
1400x2000	1530	2080	370	805	-	30°	70°	30°	240

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Alois


Model Alois

Security, practicality and esthetics. Low height bed for patients with Alzheimer's disease. Very low minimum height and maximum height up to 720 mm. Compatible with most patientlifts on the market.



Options

- 3 functions electrical bed with knee buckle
- Protective transport cover for the wooden panels
- Serum holder
- Louis Philippe wooden panels
- Remote control support
- Independent serum holder



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Eagle 620 & 625

Colour



White

The Vermeiren patientlift was designed to bring a real comfort to the patient during the transfers but also for the people who are taking care of the patient.

The patientlift Eagle 620 and Eagle 625 are equipped with the Linak HB7X engine feeding from a battery with an autonomy of approximately 45 lifts.



Optional: weighing scale

Users weight
Eagle 620 : 150 kg
Eagle 625 : 175 kg

Foldable model: easy storage




Vermeiren designed this patientlift with specific attention to functionality, safety and reliability.

Slings

SANITARY

Available in
M - L - XL (sizes in mm)


	A	B	C	D
M	780	250	600	300
L	840	260	700	360
XL	910	270	820	440



COMFORT

Available in
M - L - XL (sizes in mm)


	A	B	C	D	E
M	780	400	450	300	390
L	840	45	500	360	410
XL	910	510	560	440	420



STANDARD

Available in
M - L - XL (sizes in mm)

	A	B	C	D
M	780	400	450	300
L	840	450	500	360
XL	910	510	560	440



Quality wheels of 100 mm ensure that this patientlift can be used on almost any underground.


620
625

1740
1840

1115
1290

550
680

1035
1225

580-1600
650-1730

480
480

37
42

150
175

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Albatros

Colour



Grey

Compact patient hoist with easy usage. The Albatros model folds in a few seconds for storage and transport.



Slings

Sling sizes

Available in
M - XL

M	XL
40 - 120 kg	100 - 200 kg



High quality wheels
with seperated brakes



1040



580



940



1040 - 1730



500



44,8



150

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Alesia

Colour



Reclinable back with two controllers for the patient and the caretaker



- Solid foam seat
- Height adjustable armrests
- Retractable footrest plate
- Central brake
- Directional wheels




- Table B12 (optional)

Options



Option: Side supports



700

540

0-270

530

1330

800

12°

30°

480

800*
1330**

33

140

(*back and legrests not inclined. ** with back and legrests inclined.)

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Relaxie

Colour



Options



Headrest



Extra wide armrests movable to the side



Adjustable in height model (from 508 to 568mm)



Reclinable footrest 500 x 500 x 450 mm (LxLxH)



Table



Model Normandie

- Reclinable back 45°

Legrest will recline automatically with the back.

- Armrests adjustable in height

Model XXL

max. users weight
kg 350



Wheel version



Normandie
XXL

660
890

540
740

0-170
0-220

510
510

1200
1200

750
750

440
480

0-45°
0-40°

780
780

26
38

130
350

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Relax Provence & Languedoc


Colour



Model Languedoc

- Reclinable back 45°

Options



Retractable footrest



Table



Reclinable footrest



Headrest

Model Provence

- Fixed back

Model with toiletseat

All models

are originally equipped with height adjustable armrests.



Provence
Languedoc

660
660

540
540

0-170
0-170

510
510

1200
1200

750
750

440
440

0-45°

780
780

19
20

130
130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Ontario Liftchairs

Colour



**Warranty
2 years**

Model Ontario




The Ontario liftchair offers you optimal comfort in your own house. With a simple push on a button you can relax, watch tv or have a little nap in your Ontario liftchair.

Model Ontario 2 motors



With the Ontario 2 motors, you can separately control the backrest and the legrest.




Ontario	810	450	180	460	1120	920	500	750	52	150
Ontario II	810	450	180	460	1120	920	500	750	52,5	150

Sizes in mm, weight in kg - Deviations in sizes and weight (1,5 mm, 1,5 kg and 1,5°) are possible

Liftchairs & Montreal

Colour

Skai
taupeTissue
beigeTissue
mottled
brown

Model Vancouver



Both liftchairs can be placed very close to the wall (100 mm) and still recline completely.




Colour

Skai
blackSkai
taupe

Model Montreal



The Vancouver and Montreal model will offer you the same possibilities as the Ontario. Make your choice out of different colors and textures for an optimal match with your existing interior. The footrest inclines first making it easy to watch television without having to lay back completely.


Vancouver
Montréal
730
790450
450140
180500
5001090
1090920
920500
530650
65056
58150
150

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

286B & 4-Light

Colour



Model 286B

4-wheel foldable rollator. Equipped with a removable basket, plexi table, seat, crutch holder and anti-tangle system for the braking cables with two small hooks. Parking brakes and height adjustable handles.

Very simple and fast folding system.



800-990

700

610

590

430

11

130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Colour



Porte canne

Model 4-Light

Lightweight foldable 4-wheel aluminium rollator.

- foldable model / easy for transport
- crutch holder
- multi-position brakes
- no visible brake cables
- canvas seat
- back support belt
- storage bag



790-990

700

560

600

400

7,5

130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible



Rollator

275 & 276 & 201 Delta

Colour



C29
Grey

Model 275

Equipped with a removable basket, a seating cushion, plexiglas tray drawer, anatomic handles, seat and back. Ultra-light foldable model with progressive brakes and parking position.



275	810 - 960	760	580	590	440	8,8	120
276	1030-1150	760	580	630	440	11	120

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Model 276

Model with extra long pads to rest your elbows.



Model 201 Delta

Foldable 3-wheel narrow rollator with a storage bag, removable tablet and shopping basket. Progressive brakes and parking position.

	870-990	580	220	630	230	7	130

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Rollator 265 & 266 & 267

Colour



C27
Blue

Model 265 demountable

Tubular steel frame painted in blue.
Padded armpit supports which are adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.



Model 266 foldable

Tubular steel frame painted in blue. Padded armpit supports which are adjustable in height. 4 wheels of 75 mm of which 2 with brake.
Foldable for easy transport and storage.



Model 267 demountable

Same features as the mod.265 with the addition of the adjustment of the base height up to 150 mm (30 mm).



265	1170-1470	760	520	580	min. 300
266	1170-1470	760	520	580	min. 275
267	1115-1560	760	465-610	580	min. 300

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

268 & 269 Rollator

Colour



C27
Blue

Model 268 foldable with table

Tubular steel frame painted in blue.
Padded arm support adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.



Model 269 Foldable with lower arm support

Tubular steel frame painted in blue.
Lower arm support adjustable in height.
4 wheels of 75 mm of which 2 with brake.
Can be dismantled for easy transport and storage.

Options for walking frames

- Table (type 268)
- Lower arm supports (type 269)
- Seat



268
269

830-1130
830-1130

760
760

520
520

580
580

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Streeters


Colour



The Vermeiren Streeter was designed to ease the learning and relearning of walking. The locking system on the backwheels will prevent the Streeter from moving backwards and will encourage the user to move forward.


The system is fully reversible so that your Streeter can be used in both directions. The aluminium frame is fully collapsible which makes it excellent for transport and storage. 4 different colors represent the 4 different sizes.

Model Small



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


Model Medium



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Streeters

Model Large



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Model X-Large




options



Rubbers



Directional wheels



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Toiletoirs

Colour frame



C27
Blue

Model 139B

Commode with wheels and removable arm- and footrests.



570 420 520 420 410 12,9 100

Colour frame



C27
Blue

Model 9139



585 450 490 420 450 13,8 100

Colour frame


Model 9060



C29
Grey



570 470 500 470 360 7 120



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


Toilechairs

9062 & 9063 & 9062XXL

Colour frame



Model 9062




Colour frame




Model 9062XXL

XXL commode with a seat width of 600 mm and max. users weight of 175 kg.




Colour frame



Model 9063

Height adjustable chair.

Height adjustment:
from 410 to 560 mm



Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Serie 9300 9300

Colour frame Colour upholstery



C29
Grey



C37
Dark blue




Modello 9300 big rear wheels

Anatomical shaped back, reclinable and removable middle part of the seat, removable rear wheels 24x1 3/8 PU, directional front wheels 200x50 PU soft, manual brakes on the rear wheels.
Seat width: 430 - 460 - 500 - 530 mm

Accessories for the commode chairs

- Headrest L55
- Headrest L06
- Spreader B22
- Side supports L04
- Wooden table B12
- Plexi table B15



L06



L55



B22

Model 9301 directional rear wheels

Anatomical shaped back, reclinable and removable middle part of the seat, removable front wheels 24x1 3/8 PU, directional rear wheels 200x50 PU soft, manual brakes on the front wheels.
Seat width: 430 - 460 - 500 - 530 mm



Serie 9300 9300




Model 9302 wheels of 200x50 mm

Anatomical shaped back, reclinable and removable middle part of the seat, 4 wheels of 200x50 PU with two directionals at the rear, manual brakes on the front wheels.
Seat width: 430 - 460 - 500 - 530 mm

Model 9303 wheels of 100 mm

Anatomical shaped back, reclinable and removable middle part of the seat, 4 wheels of 100 mm with brakes on the rear wheels.
Seat width: 430 - 460 - 500 - 530 mm



9300	660	430-530	260	500	580	450	1070	1100	25,5	100
9301	660	430-530	260	500	580	450	1070	890	26,5	
9302	620	430-530	260	500	580	450	1070	840	24,5	
9303	560	430-530	260	500	580	450	1070	850	23,5	

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible





tricycles
Tricycles

Tricycles 2102 - 2103 - 2115 - 2117



Model 2102

Freedom of movement is important for everyone and almost equals independence. Furthermore, mobility and even more so cycling are important parts of rehabilitation.

That's why Vermeiren offers a wide range of tricycles. Five different models, each with its own characteristics and features, the result of more than 55 years of research and experience.



Model 2115



Model 2103



Model 2117

The table on the next page gives you a full overview of the various models, including detailed specifications and information about the available colours and options.



2102	1100	660	500 - 700	710 - 840	310
2103	1190	680	550 - 750	760 - 910	390
2115	1330	680	550 - 750	760 - 910	390
2117	1630	680	720 - 950	900-1050	480
2119	1820	810	820-1050	1030-1150	580

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

2119 Options

Model 2119



Options

- 01 Neutral
- 02 Back pedal brake
- 04 Front and rear light
- 05 Drumbrake, built into front wheel axle
- 06 Fixed sprocket
- 08 T-shaped saddle pin allowing the saddle to be adjusted back and forth
- 09 Velcro pedals
- 10 Footclasp
- 11 Velcro footplates (small model)
- 12 Velcro footplates (big model)
- 14 Broad hip support with belt
- 15 Crutch holder
- 16 Rectangular steeringwheel (small model)
- 17 Rectangular steeringwheel (big model)
- 18 Kilometer counter
- 19 Adjustable pedalbars
- 20 Footplate with legclamps
- 21 Parking brake
- 22 Back pedal brake, 3-speed
- 23 Differential
- 24 Accompanying bar
- 25 Accompanying bar with safety holder



04



04



05



08



09



10



11



14



15



17



18



19



20



24



25

Model	01	02	04	05	06	08	09	10	11	12	14	15	16	17	18	19	20	21	22	23	24	Colour
2102	<input type="checkbox"/>	<input type="checkbox"/>			s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>				<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	red and black
2103	<input type="checkbox"/>	<input type="checkbox"/>			s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>				<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	yellow and black
2115	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	yellow and black
2117	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			blue
2119	<input type="checkbox"/>	<input type="checkbox"/>	s	s	s	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	s		grey

s - Standard ☐ - Option - - impossible





easy life
Bathroom products

Bathing seats

Colour



Blue White

Easy life

All our *Easy life* products are ideal for usage in the shower, bathtub, bathroom or other humid environments.

The specially designed rubbers at the end of each leg will create a safe stability for the bench and avoid the sliding of the bench.

Jill



Grace

Removable shower seat for the wall. Easy installation.

Jinny



Sidney

Bathroom seat with toilet opening

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible


Bathing seats



Lilly

Easy life

The transfer benches Kate and Katy will facilitate the transfer from a chair to your bath. Like all the products in the range, they are equipped with a non-skid surface. The handle mounts right or left of the bench depending on your needs.



Kate



NEW

Paris

Colour




Katy

Transfer bench with bucket



Colour



Melbourne

Bath chair with toilet opening

Sizes in mm, weight in kg - Deviations in sizes and weight (1,5 mm, 1,5 kg and 1,5°) are possible

Toilets Easy life

Colour



White

Stacy

Commode
Raised toilet seat
Toilet frame
Armrests can be lowered for easy transfer.



Tina

The Vermeiren raised toilet seat Tina was designed to raise the height of a toilet to a comfortable level, making it easier to sit down or get up.

width: 520 mm
depth: 430 mm
height: 380 mm
seat width: 450 mm
adjustable armrests from 5 to 90 mm
height without armrests: 110 mm
weight: 2,5 kg
users weight: 136 kg



Liddy

Toilet frame adjustable in height and ease of use.
width : 560 mm
depth : 470 mm
height : 630 - 740 mm
weight : 1,45 kg
users weight : 113 kg

Installation without drilling.



General life

Colour




Ivy



Functional and sturdy bathtub rail providing extra support and safety whilst entering or exiting your bathtub.

Easy to install and applicable on most types of bathtubs due to its adjustable width (from 60 to 115 mm).

Iris




Izzie

The bathroom bars are available in different lengths from 300 to 800 mm. With exception of the Izzie R and Izzie L, which are bars in L-shape. They are designed as general support elements. They are useable everywhere, including humid environments.

Izzie 300	400	450	600	800
0,4	0,45	0,50	0,55	0,75 kg

Izzie R




Sizes in mm, weight in kg - Deviations in sizes and weight (1,5 mm, 1,5 kg and 1,5°) are possible

General life

Foldable canes



Ashley Tall & Medium (adjustable height)



Easy life

Vermeiren offers a wide variety in canes for in- and outdoor usage. From stylish foldable models to fixed models.

Find the color and design that suits you most and have more enjoyment of your daily walks.



Zoë


Fixed canes



Nicky



The helping hand Nicky is designed for daily assistance to grab or pick up small objects on the floor. The front of the helping hand Nicky is equipped with a small magnet. It can help you to pick up small magnetic objects such as keys, coins or cans.



thickness: 12 mm
length: 760 mm
height: 110 mm
weight: 0,15 kg

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Walkers Easy life

Lyn (fixed foldable model)

width: 580 mm
depth: 480 mm
folded: 110 mm
height: 810 - 990 mm
weight: 2,8 kg
users weight: 113 kg



Lyna (articulating foldable model)

width: 580 mm
depth: 480 mm
folded: 110 mm
height: 810 - 990 mm
weight: 2,8 kg
users weight: 113 kg



Lindsey (fixed foldable model)

Model "sit/stand-up"
width: 630 mm
depth: 540 mm
folded: 120 mm
height: 790 - 970 mm
weight: 2,75 kg
users weight: 113 kg



convenient height above the toilet

Sizes in mm, weight in kg - Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible

Compatibility table

	Sagitta	28	708	708 kids	Jazz	Jazz 30°	Inovys	Serenys	917	925	Springer	Express	Sunny	Tracer kids	Tracer (GT)	Rapido	Squod	Squod S.U.	Ascend	Forest	Eclips	Eclips 30°	Eclips X4	Eclips X4 kids
B02	-	•	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-	•	-
B03	-	•	•	•	•	•	-	•	-	-	-	•	-	-	-	•	-	-	-	-	•	•	•	-
B05	-	•	•	-	•	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•	•
B06	•	•	•	•	•	•	-	•	-	•	•	•	•	•	•	•	•	-	•	•	•	•	•	•
BA6	•	•	•	•	•	•	-	-	-	-	•	•	-	•	•	•	•	-	•	•	•	•	•	•
BK6	•	•	•	•	•	•	-	•	-	-	•	•	•	•	•	•	•	•	•	•	•	•	•	•
BL6	-	•	•	•	•	•	-	-	-	-	-	•	-	•	•	•	•	-	•	•	•	•	•	•
BL7	-	•	•	•	•	•	-	-	-	-	-	-	-	-	-	•	•	-	•	•	•	•	•	-
BK7	-	•	•	•	•	•	-	-	-	-	-	-	-	-	-	•	•	-	•	•	•	•	•	•
BN6	-	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B07	-	•	•	•	-	-	-	-	-	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-
BZ7	-	•	•	•	•	•	•	•	-	-	-	•	-	-	•	•	•	-	•	•	•	•	•	-
B08	-	•	•	•	-	-	-	-	-	•	-	•	•	-	-	-	-	-	-	-	-	-	-	-
BZ8	-	•	•	•	•	•	•	-	•	-	•	•	-	-	•	•	•	-	•	•	•	•	•	•
B09	-	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
BZ9	-	•	•	•	-	•	-	-	•	-	-	-	-	-	-	-	-	-	-	-	•	•	•	-
B10	-	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
BZ10	-	•	•	•	-	•	-	-	•	-	-	-	-	-	-	-	-	-	-	-	•	•	•	-
B12	-	•	•	•	•	•	•	•	•	-	-	•	•	-	•	•	•	-	•	-	•	•	•	•
B13	-	•	•	•	•	•	•	•	•	-	•	•	•	•	•	•	•	-	•	•	•	•	•	•
B14	-	•	•	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•	-
B15	-	•	•	-	•	•	-	•	-	-	•	-	-	•	-	-	-	-	-	•	•	-	•	-
B17	-	•	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B18	-	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B19	-	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B20	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B21	-	•	•	•	•	•	•	•	-	•	•	•	•	•	•	•	•	-	•	•	•	•	•	•
B22	-	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	-	•	•	•	-	-	-	-
B23	•	•	•	•	•	•	•	•	-	•	•	•	•	•	•	•	•	-	•	•	•	•	•	•
B24	-	•	•	•	•	•	•	•	-	•	-	•	•	•	•	•	•	-	•	•	•	•	•	•
B25	-	•	•	•	•	•	•	•	-	•	-	•	•	•	•	•	•	-	•	•	•	•	•	•
B26	-	•	•	•	•	•	•	•	-	•	-	•	•	•	•	•	•	-	•	•	•	•	•	•
B27	-	•	•	•	•	•	•	•	-	•	-	•	•	•	•	•	•	-	•	•	•	•	•	•
B28	-	•	•	•	•	•	-	-	-	•	•	-	•	•	•	-	-	-	-	•	•	-	•	-
B30	-	•	•	•	•	•	•	•	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•	•
B31	-	•	•	•	•	•	-	•	-	-	•	•	-	•	•	•	•	•	•	•	•	•	•	-
B32	-	•	•	•	•	•	•	•	•	•	-	•	•	-	•	•	•	-	•	-	•	•	•	-
B37	-	•	•	-	•	•	-	•	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•	-
B40	-	•	•	•	•	•	•	•	•	•	-	-	•	-	•	•	•	-	-	-	•	•	•	•
B41	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B42	•	•	•	•	•	•	•	•	-	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
B45	-	•	•	•	•	•	-	-	-	•	-	•	•	-	•	•	•	-	-	-	•	•	•	•
B46	-	•	•	•	•	•	-	-	-	•	-	•	•	-	•	•	•	-	-	-	•	•	•	•
	Sagitta	28	708	708 kids	Jazz	Jazz 30°	Inovys	Serenys	917	925	Springer	Express	Sunny	Tracer kids	Tracer (GT)	Rapido	Squod	Squod S.U.	Ascend	Forest	Eclips	Eclips 30°	Eclips X4	Eclips X4 kids

Compatibility table

	Sagitta	28	708	708 kids	Jazz	Jazz 30°	Innovys	Serenys	917	925	Springer	Express	Sunny	Tracer kids	Tracer (GT)	Rapido	Squod	Squod S.U.	Ascend	Forest	Eclips	Eclips 30°	Eclips X4	Eclips X4 kids
B47	-	●	●	-	●	●	-	-	-	-	-	●	●	-	-	●	●	-	-	-	●	●	●	-
B52	-	●	●	●	●	●	●	●	●	-	-	-	-	-	●	-	-	-	-	-	●	-	●	-
B57	●	●	●	●	●	●	●	●	-	●	●	●	●	●	●	●	●	-	●	●	●	●	●	●
B58	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
B59	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	●	●	●	●	●	-
B63	-	●	●	-	●	●	●	●	●	●	-	●	●	-	●	●	●	-	●	-	●	●	●	●
B64	-	●	●	-	●	●	●	●	●	●	-	●	●	-	●	●	●	-	●	-	●	●	●	●
B66	-	●	●	-	●	●	●	-	-	-	-	●	-	-	-	-	-	-	-	-	●	●	●	●
B69	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	-
B73	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B74	-	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●
B78	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
B79	-	●	●	-	●	●	-	-	-	●	-	●	-	-	-	-	-	-	-	-	●	●	●	-
B80	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●
B82	●	-	-	-	-	-	-	-	-	●	●	-	-	●	-	-	-	●	-	●	-	-	●	●
B83	-	●	●	-	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	●
B84	-	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●
B85	●	●	●	●	●	●	●	●	-	●	-	-	-	-	-	-	-	-	-	-	●	●	●	●
B90	-	●	●	-	●	●	●	-	-	●	-	●	●	-	-	●	●	-	-	-	●	●	●	●
B91	-	●	●	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	-
B92	-	●	●	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	-
B94	-	●	●	-	●	●	●	●	●	●	-	●	●	-	●	●	●	-	●	-	●	●	●	●
B99	-	●	●	-	●	●	-	-	-	-	-	●	-	-	●	●	●	-	-	●	●	●	●	-
B101	-	●	●	●	●	●	-	●	-	-	●	●	●	●	●	●	●	●	-	●	●	●	●	-
B102	-	●	●	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B112	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B115	-	●	●	-	●	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	-
AMP	-	●	●	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	-
S1	-	●	●	●	●	●	●	-	●	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●
HEM2	-	●	●	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	-
T30	-	●	●	-	●	●	-	●	●	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●
L04	-	●	●	-	●	●	●	●	-	●	-	-	●	●	●	●	●	-	●	●	●	●	●	●
L07	-	●	●	-	●	●	●	●	-	●	-	-	-	-	-	-	-	-	-	-	●	●	●	●
L14	-	●	●	-	●	●	-	-	-	●	-	●	●	-	-	●	●	-	-	-	●	●	●	●
L15	-	●	●	-	●	●	-	-	-	●	-	●	●	-	-	●	●	-	-	-	●	●	●	●
L24	-	●	●	-	●	●	-	-	-	-	-	●	-	-	●	-	-	-	-	-	●	●	●	●
L25	-	●	●	-	●	●	-	-	-	-	-	●	-	-	●	-	-	-	-	-	●	●	●	●
L55	-	●	●	-	●	●	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Sagitta	28	708	708 kids	Jazz	Jazz 30°	Innovys	Serenys	917	925	Springer	Express	Sunny	Tracer kids	Tracer (GT)	Rapido	Squod	Squod S.U.	Ascend	Forest	Eclips	Eclips 30°	Eclips X4	Eclips X4 kids

Subject to technical alterations

● possible - not possible

Notes

Handwriting practice lines consisting of 20 horizontal dotted lines.

All rights reserved, including translation.

No part of this manual may be reproduced in any form whatsoever (print, photocopy, microfilm or any other process) without written permission of the publisher, or processed, duplicated, or distributed by using electronic systems.

© N.V. Vermeiren N.V. 2012

VERMEIREN N.V.

Vermeirenplein 1-15

B2920 - Kalmthout

Belgium

Tel.: +32 (0) 3 620 20 20

Fax: +32 (0) 0 666 48 94

Mail: info@vermeiren.be



For more information, contact your local distributor.

